




ENT - Nasal

“Empowering Precision in Every Incision.”





Thomson surgicals instruments possesses several key features that contribute to its effectiveness in various medical procedures:

Precision and Accuracy:

Instruments designed for precise and accurate performance to enhance surgical outcomes.

Ergonomic Design:

Ergonomically designed handles for comfortable grip and reduced hand fatigue during long surgeries.

Durability and Material Quality:

Use of high-quality materials such as stainless steel grade 420, 410, 316, and 304 or titanium medical grade, Copper, Brass and Hylam Material for durability and resistance to corrosion.

Customization:

Offering customizable options based on surgeon preferences or specific surgical procedures.

Innovative Coatings:

Integration of advanced coatings for reduced friction, enhanced visibility, long life and ease of cleaning.

Minimally Invasive Solutions:

Developing instruments suitable for minimally invasive procedures, promoting smaller incisions and faster recovery times.


Modularity:

Modular instrument systems that allow for flexibility and adaptability in various surgical settings.

Easy Maintenance and Sterilization:

Designing instruments that are easy to clean, sterilize, and maintain to ensure optimal performance and infection control.

The combination of these features ensures that the Thomson surgicals instruments contribute positively to the surgical process, providing optimal exposure, accessibility, and safety during procedures.





CINELLI Lower Lateral Scissors

Product Number	Balde	Length (in)
EN0001	STRONG Cd	4 3/8

COTTLE Bulldog Scissors

Product Number	Tip	Blade	Length (in)
EN0002	B/B	Cd	4 1/2





COTTLE-KNAPP Plastic Scissors

Product Number	Tip	Blade	Length (in)
EN0003	B/B	ST	4 1/8
EN0004	B/B	Cd	4 1/8

COTTLE-MASING Plastic Scissors

Product Number	Tip	Blade	Length (in)
EN0005	S/S	Cd	4 1/8
EN0006	B/B	Cd	4 1/8





EN0007


DAVIS Rhytidectomy Scissors

Product Number	Blade	Length (in)
EN0007	Cd	6 3/4

EN0007

Rhinoplasty Scissors

Product Number	Tip	Blade	Length (in)
EN0008	B/B	Flat, Cd	5 1/4




EN0008

EN0008

Nasal Scissors, Pediatric

Product Number	Blade	Blade Size (mm)	Length (in)
EN0009	ST	6	4
EN0010	Cd-RIGHT	5.5	4
EN0011	Cd-LEFT	5.5	4



EN0009 - EN0011



EN0009




EN0010



EN0011

Nasal Scissors, Adult

Product Number	Blade	Blade Size (mm)	Working Length (in)
EN0012	ST	11	5 1/8
EN0013	RIGHT	11	5 1/8
EN0014	LEFT	11	5 1/8



EN0012 - EN0014



EN0012



EN0013



EN0014

LAKESIDE Nasal Scissors


Product Number	Blade	Length (in)
EN0015	ST	7 1/4



EN0015

HEYMANN Nasal Scissors

Product Number	Tip	Blade	Length (in)
EN0016	B/B	ST	7 1/2



EN0016

EN0016

HEYMANN-KNIGHT Nasal Scissors

Product Number	Tip	Blade	Length (in)
EN0017	B/B	ST	6 7/8




EN0017

EN0017


KNIGHT Nasal Scissors, Lightweight, Delicate

Product Number	Tip	Blade	Length (in)
EN0018	DEL	ST	6 3/8



JACKSON Nasal Bone Scissors


Product Number	Blade	Length (in)
EN0019	ANG	8 1/2



Dorsal Scissors


Product Number	Tip	Blade	Length (in)
EN0020	BLUNT	ST	4 3/4





EN0021

EN0021



EN0022


EN0022

FOMON Dorsal Scissors

Product Number	Tip	Blade	Handle	Length (in)
EN0021	B/B	SERR	ANG	5 3/4
EN0022	B/B	SERR	ANG	5 3/4

FOMON Lower Lateral Scissors

Product Number	Tip	Blade	Length (in)
EN0023	B/B	HEAVY-Cd	4 1/2
EN0024	B/B	HEAVY-Cd	5 1/4



EN0023

EN0024

EN0023 - EN004



EN0025


EN0025

FOMON Saber-Back Scissors

Product Number	Tip	Blade	Length (in)
EN0025	B/B	Cd, SERR	5


GOLDMAN Septum Scissors

Product Number	Blade	Length (in)
EN0026	ST, SERR	5 3/4



COTTLE Septum Scissors


Product Number	Blade	Length (in)
EN0027	HEAVY, ANG	7 1/2





GORNEY Turbinate Scissors, Double Action

Product Number	Jaw	Length (in)
EN0028	SERR	8



EN0028


EN0028

Becker Septum Scissors, Double Joint


Product number	Blade	Length (in)
EN0029	SERR	7



EN0029



EN0029



EN0030


EN0030

MARTIN Cartilage Scissors

Product Number	Blade	Length (in)
EN0030	Cd, SERR	8

METZENBAUM-LIPSETT Scissors

Product Number	Tip	Blade	Length (in)
EN0031	B/B	Cd	5 1/2




EN0031

EN0031

PECK-JOSEPH Scissors

Product Number	Blade	Tip Size (mm)	Length (in)
EN0032	ST, SHARP	42	5 3/4




EN0032

EN0032

KAYE Facelift Scissors



EN0033




EN0033

Product Number	Tip	Blade	Length (in)
EN0033	SEMI-SHARP	Cd, SERR	6
EN0034	SEMI-SHARP	Cd, SERR	7 1/2



EN0033



EN0034

KILNER Scissors

Product Number	Tip	Blade	Length (in)
EN0035	FLAT	Cd	5
EN0036	FLAT	Cd	6



EN0035 - EN0036




EN0035



EN0036

JAMESON Tenotomy Scissors, Delicate

Product Number	Tip	Blade	Length (in)
EN0037	B/B	Cd	6




EN0037

EN0037

REYNOLDS-JAMESON Tenotomy Scissors, Delicate

Product Number	Tip	Blade	Handle	Length (in)
EN0038	B/B	Cd	SQUARE	5 1/2




EN0038

EN0038

REYNOLDS Tenotomy Scissors, Delicate

Product Number	Tip	Blade	Length (in)
EN0039	B/B	Cd	5 3/4
EN0040	B/B	Cd	6 7/8




EN0039

EN0039 - EN0040



EN0040

Iris Scissors, Delicate



Product Number	Tip	Blade	Handle	Length (in)
EN0041	S/S	ST	SQUARE	4 1/2
EN0042	S/S	Cd	SQUARE	4 1/2

EN0041

EN0041 - EN0042



EN0042

STEVENS Tenotomy Scissors, Delicate

Product Number	Tip	Blade	Handle	Length (in)
EN0043	S/S	ST	RIBBON	4
EN0044	S/S	Cd	RIBBON	4
EN0045	B/B	Cd	RIBBON	4 1/4




EN0043

EN0043 - EN0044



EN0044



EN0045

EN0045


STEVENS Tenotomy Scissors, Delicate




EN0046 - EN0049

Product Number	Tip	Blade	Handle	Length (in)
EN0046	S/S	ST	SQUARE	4 1/2
EN0047	B/B	ST	SQUARE	4 1/2
EN0048	S/S	Cd	SQUARE	4 1/2
EN0049	B/B	Cd	SQUARE	4 1/2

EN0046



EN0047



EN0048



EN0049

JOSEPH Scissors, Delicate

Product Number	Tip	Blade	Handle	Length (in)
EN0050	S/S	ST, SERR	ROUND	5 1/2
EN0051	S/S	Cd, SERR	ROUND	5 1/2
EN0052	S/S	ST	ROUND	6
EN0053	S/S	Cd	ROUND	6
EN0054	S/S	ST	SQUARE	6
EN0055	S/S	Cd	SQUARE	6



EN0050 - EN0051



EN0050




EN0051



EN0052 - EN0053




EN0052



EN0053



EN0054 - EN0055




EN0054




EN0055

STEVENS Tenotomy Scissors

Product Number	Tip	Blade	Handle	Length (in)
EN0056	S/S	ST	ROUND	4 1/4
EN0057	B/B	ST	ROUND	4 1/4
EN0058	B/B	ST	ROUND	5
EN0059	B/B	ST	ROUND	5 3/4
EN0060	S/S	Cd	ROUND	4 1/4
EN0061	B/B	Cd	ROUND	4 1/4
EN0062	B/B	Cd	ROUND	5
EN0063	B/B	Cd	ROUND	5 3/4



EN0056 - EN0059



EN0056




EN0057




EN0058




EN0059




EN0060 - EN0063



EN0060



EN0061




EN0062



EN0063


STEVENS Tenotomy Scissors

Product Number	Tip	Blade	Handle	Length (in)
EN0064	B/B	ST	RIBBON LIGHT	3 3/4
EN0065	B/B	ST	RIBBON HEAVY	3 7/8
EN0066	B/B	Cd	RIBBON HEAVY	3 3/4




EN0064

EN0064



EN0065

EN0065



EN0066

EN0066

Strabismus Scissors, Delicate

Product Number	Tip	Blade	Handle	Length (in)
EN0067	B/B	ST	SQUARE	4 1/2
EN0068	B/B	Cd	SQUARE	4 1/2



EN0067

EN0068

EN0067 - EN0068



FOMON Dorsal Scissors with tungsten carbide edges

Product Number	Tip	Blade Angle	Length (in)
EN0069	B/B	30°	5 1/4




EN0069

EN0069

Iris Scissors, Delicate with tungsten carbide edges


Product Number	Tip	Blade	Length (in)
EN0070	S/S	ST	4 1/4
EN0071	S/S	Cd	4 1/4



EN0070




EN0071



EN0070 - EN0071

Iris Scissors, Wavecut with tungsten carbide edges

Product Number	Tip	Blade	Length (in)
EN0072	S/S	Cd	4 1/4



EN0072

EN0072

Strabismus Scissors, Delicate with tungsten carbide edges



EN0073 - EN0074

Product Number	Tip	Blade	Length (in)
EN0073	B/B	ST	4 1/2
EN0074	B/B	Cd	4 1/2




EN0073



EN0074

GOLDMAN-FOX Tenotomy Scissors, Delicate with tungsten carbide edges

Product Number	Tip	Blade	Length (in)
EN0075	S/S	CD, 1/SERR	5 1/8




EN0075


EN0075

STEVENS SuperCut Scissors

Product Number	Tip	Tip Size (mm)	Blade	Length (in)
GS0228	B/B	12	ST	4 1/8
GS0229	B/B	12	Cd	4 1/8
GS0230	B/B	13	ST	4 1/4
GS0231	B/B	13	Cd	4 1/4



GS0228



GS0230



GS0231



GS0229

FOMON SuperCut Scissors

Product Number	Tip	Blade	Length (in)
EN0076	B/B	Cd	4 1/2



EN0076

EN0076




EN0077

Septum Needle W/ Luer Lock, Attachment

Product Number	Tip Size (mm)	Length (in)
EN0077	DIAM 0.8	3 1/2

Antrum Trocar Needle

Product Number	Tip	Size	Length (in)
EN0078	BEVELED STYLET	14 GA	4 1/4
EN0079	BEVELED STYLET	16 GA	4 1/4
EN0080	BEVELED STYLET	18 GA	4 1/4



EN0078 - EN0080



EN0078 - EN0080

Septum Needle

Product Number	Tip	Size	Length (in)
EN0081	ANG	22	5 3/8



EN0081

COAKLEY Antrum Trocar, with Irrigation Extension

Product Number	Length (in)
EN0082	6 1/4



EN0082

EN0082

Sinoscopy Trocar and Cannula



product Number	Tip Size (mm)	Length (in)
EN0083	4	6 1/2
EN0084	5	6 1/2


OCHSNER Trocar

Product Number	Tip size(FR)	Length(in)
GS1856	10	6
GS1857	12	6
GS1858	14	6
GS1859	16	6
GS1860	18	6
GS1861	20	6
GS1862	22	6
GS1863	24	6
GS1864	26	6
GS1865	27	6
GS1866	28	6



BROWN Needle Holder, Heavy

Product Number	Tip Size (mm)	Jaw	Length (in)
EN0085	13	SERR	5 1/4



EN0085

EN0085

SOFFERMAN Needle Holder

Product Number	Tip Size (mm)	Jaw	Working Length (in)
EN0086	5	ST, SERR	3 5/8




EN0086

EN0086



GS0482



GS0482

CALLISON-ADSON Tissue Forceps, W/ Platform


Product Number	Teeth	Length (in)
GS0482	1 x 2	5

LUCAE Bayonet Forceps

Product Number	Tip	Teeth	Length (in)
GS0438	SERR	-	5
GS0439	SMOOTH	1 x 2	5
EN0087	SERR	-	6 3/8
EN0088	DEL, SERR	-	5 3/4



GS0438 - GS0439



GS0438



GS0439



EN0087



EN0087



EN0088




EN0088

JANSEN Bayonet Forceps

Product Number	Tip	Length (in)
GS0445	SERR	6
GS0446	SERR	6



GS0445 - GS0446



GS0445



GS0446

GRUENWALD Bayonet Forceps




Product Number	Tip	Length (in)
GS0440	SERR	8

JACOBSON Forceps

Product Number	Tip	Length (in)
EN0089	SMOOTH, SHARP	6 1/4



Bayonet Forceps




Product Number	Tip	Teeth	Length (in)
GS0447	SMOOTH	-	7
GS0448	SMOOTH	1 x 2	8

BROWN-BAHNSON Bayonet Forceps

Product Number	Tip	Length (in)
GS0444	SERR	6



ADSON Bayonet Forceps




Product Number	Tip	Length (in)
EN0090	SERR	6 7/8




WILDE Forceps

Product Number	Tip	Length (in)
GS0457	ANG/SERR	4



TROELTSCH Bayonet Forceps

Product Number	Tip	Length (in)
GS0454	ANG/SERR	4
GS0455	ANG/SERR	5
GS0456	ANG/SERR	7




HARTMANN Alligator Forceps

Product Number	Tip	Length (in)
EN0091	ANG/SERR	7 1/2



HARTMANN Nasal Dressing Forceps


Product Number	Length (in)
EN0092	6 7/8
EN0093	7
EN0094	7 7/8
EN0095	8 1/2





WEST Nasal Dressing Forceps

Product Number	Tip	Length (in)
EN0096	SERR	7 7/8



COTTLE Columella Clamp, W/ Fixation Screw

Product Number	Length (in)
EN0097	4



FERRIS-SMITH Fragment Forceps

Product Number	Tip Size (mm)	Working Length (in)	Length (in)
EN0098	8.2 x 6	4 1/4	7 7/8



OSTRUM Antrum Backbiting Punch Forceps

Product Number	Jaw Size (mm)	Jaw	Working Length (in)	Length (in)
EN0099	3 x 7	ST	4 1/4	7 3/4
EN0100	3 x 7	ANG-LEFT	4 1/4	7 3/4
EN0101	3 x 7	ANG-RIGHT	4 1/4	7 3/4



EN0099



EN0100



EN0101



EN0099 - EN0101

BRUENINGS Septum Forceps

Product Number	Size	Length (in)
EN0102	SMALL	7 1/2
EN0103	MEDIUM	7 1/2
EN0104	LARGE	7 1/2
EN0105	X-LARGE	7 1/2



EN0102 - EN0105



EN0102



EN0104



EN0103



EN0105

HEYMANN-KNIGHT Septum Forceps

Product Number	Jaw	Length (in)
EN0106	OVAL	7 1/2



EN0106



EN0106




JANSEN-MIDDLETON Septum Forceps

Product Number	Jaw	Length (in)
EN0107	SPOON	7 7/8
EN0108	THRU-CUT	8 1/4



EN0107



EN0107



EN0108



EN0108

GRUENWALD Nasal Forceps

Product Number	Tip Size (mm)	Tip Size (#)	Working Length (in)	Length (in)
EN0109	2.5 x 9	1	4	7 1/4
EN0110	3 x 10	2	4	7 1/4
EN0111	3.4 x 10	3	4	7 1/4
EN0112	3.5 x 10	3	4	7 1/4



EN0109 - EN0112

EN0109

EN0111

EN0110

EN0112


GRUENWALD-HENKE Nasal Thru-Cut Rongeur

Product Number	Jaw	Working Length (in)
EN0113	SMALL	4 3/4
EN0114	LARGE	4 3/4




STAMMBERGER Circular Punch Forceps

Product Number	Tip	Tip Size (mm)	Working Length (in)	Length (in)
EN0115	ST	3.5	7 1/2	9 3/4
EN0116	ST	4.5	7 1/2	9 3/4
EN0117	70° UP	3.5	7 1/2	9 5/8
EN0118	70° UP	4.5	7 1/2	10 1/4



Giraffe Frontal Sinus Forceps

Product Number	Jaw	Cup Size (mm)	Working Length (in)	Length (in)
EN0119	ANG 45° VERTICAL CUTTING	2	5 1/2	8 7/8
EN0120	ANG 45° HORIZONTAL CUTTING	2	5 1/2	8 7/8
EN0121	ANG 90° VERTICAL CUTTING	2	5 1/8	8 3/8
EN0122	ANG 90° HORIZONTAL CUTTING	3	5 1/8	8 3/8
EN0123	ANG 45° VERTICAL CUTTING	3	5 1/2	8 7/8
EN0124	ANG 70° VERTICAL CUTTING DOUBLE ACTION	3	5 1/8	8 1/4
EN0125	ANG 90° HORIZONTAL CUTTING	3	5 1/2	8 7/8
EN0126	ANG 70° HORIZONTAL CUTTING DOUBLE ACTION	4	5 1/8	8 1/4
EN0127	ANG 100° HORIZONTAL CUTTING DOUBLE ACTION	4	4 5/8	8 1/4



Giraffe Thru-Cut Sinus Forceps


Product Number	Jaw	Cup Size (mm)	Working Length (in)
EN0128	ANG 70° VERTICAL CUTTING	3 x 6	4 3/4






Antrum Punch Forceps

Product Number	Jaw	Tip	Cup Size (mm)	Working Length (in)	Length (in)
EN0129	ANG 45° BACKWARD 3MM	3	4 x 5	5 1/2	8 1/4
EN0130	ANG 45° FORWARD 3MM	3	4 x 5	5 1/2	8 1/4




MCKENTY Sphenoid Punch Forceps, Rotatable

Product Number	Jaw	Working Length (in)	Length (in)
EN0131	1.5 x 4	5 3/4	7




KERRISON Nasal Rongeur

Product Number	Jaw Size (mm)	Jaw	Working Length (in)
EN0132	4.0	90° UP	3 1/4
EN0133	5.0	90° UP	3 3/8



HAJEK-KOFLER Punch


Product Number	Jaw Size (mm)	Jaw	Working Length (in)
EN0134	3.3	90° UP	5 1/2
EN0135	3.4	90° DOWN	5 1/2
EN0136	4.0	90° UP	5 1/2






HAJEK-KOFLER Rotatable Punch, (continued)

Product Number	Jaw Size (mm)	Jaw	Working Length (in)
EN0137	4.0	90° UP	6 3/4




KOFLER Sinus Punch, Rotatable

Product Number	Jaw Size (mm)	Jaw	Working Length (in)
EN0138	3.0	90° UP	5 3/4



TAKAHASHI Nasal Forceps

Product Number	Jaw Size (mm)	Jaw	Working Length (in)
EN0139	2.5	ST	4 3/8
EN0140	3.0	ST	4 3/4
EN0141	3.0	ST	6 3/8
EN0142	4.0	ST	4 3/8
EN0143	4.0	ST	6 3/8
EN0144	3.0	45° UP	4 3/4



LOVE-GRUENWALD Rongeur


Product Number	Jaw Size (mm)	Jaw	Working Length (in)
EN0145	2 x 10	ST	5 1/4
EN0146	3 x 10	ST	5 1/4





SPURLING Rongeur

Product Number	Jaw Size (mm)	Jaw	Working Length (in)
EN0147	4 x 10	ST	5 1/8



EN0147



EN0147

MELTZER Adenoid Punch Forceps

Product Number	Tip Size (mm)	Size (#)	Jaw	Working Length (in)	Length (in)
EN0148	5MM BOTTOM BASKET, 4MM UPPER JAW	0	TRIANGULAR	4 1/2	7 1/4
EN0149	7 MM BOTTOM BASKET, 3 MM UPPER JAW	1	TRIANGULAR	4 1/2	7 1/4
EN0150	9MM BOTTOM BASKET, 8 MM UPPER JAW	2	TRIANGULAR	4 1/2	7 1/4
EN0151	6MM BOTTOM BASKET, 5MM UPPER BITE	1	ROUND	4 1/2	7 1/4

EN0148



EN0148 - EN0150

EN0149



EN0150



EN0151




EN0151



Nasal Biopsy Forceps

Product Number	Cup Size (mm)	Tip	Working Length (in)	Length (in)
EN0152	2.5 x 4	ST	5	7
EN0153	2.5 x 5	70° UP	4	7
EN0154	2.5 x 6	70° UP	3 1/2	7 1/4

EN0152



EN0152

EN0153



EN0153

EN0154



EN0154



Nasal Suction Forceps

Product Number	Jaw Length (mm)	Tip	Tip Size (#)	Working Length (in)	Length (in)
EN0155	2.5	ST	00	4	7 1/4
EN0156	3.5	ST	0	4	7 1/4
EN0157	4.5	ST	1	4	7 1/4
EN0158	2.5	45° ANG-UP	00	4 3/4	7 1/4
EN0159	3.5	45° ANG-UP	0	4 1/2	7 3/4



EN0155



EN0156



EN0157



EN0155 - EN0157



EN0158



EN0158




EN0159



EN0159

Pediatric Ostrum Antrum Punch Forceps


Product Number	Jaw Length (mm)	Jaw	Working Length (in)	Length (in)
EN0160	1.5 x 4	ST BACK BITE	4 1/2	8
EN0161	1.5 x 4	RIGHT SIDE BITE	4 1/2	8
EN0162	1.5 x 4	LEFT SIDE BITE	4 1/2	8





RONIS Adenoid Nasal Forceps

Product Number	Tip Size (mm)	Tip Size (#)	Jaw	Working Length (in)	Length (in)
EN0163	7 MM BASKET	1	TRIANGULAR	5 3/8	8 1/8
EN0164	8 MM BASKET	2	TRIANGULAR	5 3/8	8 1/8




STAMMBERGER Sidecutting Antrum Punch Forceps

Product Number	Jaw	Length (in)
EN0165	LEFT CUT	8 5/8
EN0166	RIGHT CUT	8 5/8



Ethmoid Forceps, Pediatric

Product Number	Jaw Length (mm)	Jaw	Working Length (in)	Length (in)
EN0167	3 x 5	ST	4	6 1/8
EN0168	3 x 5	45° ANG UP	4	6 1/8
EN0169	3 x 5	90° ANG UP	3 3/4	6 1/8



EN0167 - EN0169



EN0167



EN0168




EN0169

Rotating Ostrum Antrum Punch Forceps, Backbiting

Product Number	Jaw Size (mm)	Size	Jaw	Working Length (in)	Length (in)
EN0170	1.5 x 4	Pediatric	360° ROTATING	4 1/4	8 1/4
EN0171	3 x 7	Adult	360° ROTATING	4 1/4	8 1/4

110 mm, 4 1/4



EN0170



EN0171

EN0170 - EN0171



Sinus Cup Forceps

Product Number	Cup Size (mm)	Working Length (in)	Length (in)
EN0172	3	4 1/2	7



EN0172

EN0172

COTTLE Biting Forceps

Product Number	Jaw Length (mm)	Jaw	Length (in)
EN0173	30	ANG	7




EN0173



EN0173

ROWLAND Bone Cutting Forceps, Double Action

Product Number	Jaw Length (mm)	Length (in)
EN0174	25	7 1/4



EN0174


STAMM Bone Cutting Ronguer

Product Number	Length (in)
EN0175	6



STUEMPEL-VOSS Nasal Forceps

Product Number	Cup (mm)	Tip	Working Length (in)	Length (in)
EN0176	3 x 5	ANG UP	4 1/4	6 7/8





WEIL Nasal Forceps


Product Number	Cup (mm)	Tip	Working Length (in)	Length (in)
EN0177	5 x 8	ST	4	7 1/8
EN0178	5 x 8	ANG UP	4	7 1/8




WEIL-BLAKESLEY Nasal Forceps

Product Number	Jaw Size (mm)	Tip	Tip Size (#)	Working Length (in)
EN0179	3.0	ST	0	4 3/4
EN0180	3.5	ST	0	4
EN0181	3.6	ST	0	4 3/4
EN0182	4	ST	1	4 1/2
EN0183	4.2	ST	1	4 3/4
EN0184	4.8	ST	2	4 3/4
EN0185	5	ST	2	4 1/2
EN0186	5.6	ST	3	4 3/4
EN0187	3.6	ANG 45° UP	0	4 1/2
EN0188	4	ANG 45° UP	1	4 1/8
EN0189	4.2	ANG 45° UP	1	4 1/2
EN0190	4.8	ANG 45° UP	2	4 1/2
EN0191	4	ANG 90° UP	1	4

120 mm, 4 1/4



EN0179 - EN0186



110 mm, 4 1/4



EN00187 - EN0190



EN0191




EN0191




BLAKESLEY Thru-Cut Forceps

Product Number	Size	Jaw Size (mm)	Tip	Jaw Size (#)	Working Length (in)	Length (in)
EN0192	PEDIATRIC	2.5	ST	00	4 1/4	7
EN0193	PEDIATRIC	2.5	ANG	00	4 1/4	7
EN0194	-	3.5	ST	0	4 1/2	7 3/4
EN0195	-	2.5	ANG UP	00	4 3/8	7 1/4




BLAKESLEY-WILDE Nasal Forceps

Product Number	Jaw Length (mm)	Tip	Length (in)
EN0196	10 x 3.5	ANG 90° UP	4 1/4
EN0197	13.5 x 3.8	ANG 90° UP	4 1/4
EN0198	14. 3 x 5	ANG 90° UP	4 1/4



WILDE Ethmoid Nasal Forceps

Product Number	Jaw Length (mm)	Tip	Jaw Size (#)	Working Length (in)
EN0199	3.5 x 12	ANG 45° UP	0	4 3/8





WATSON-WILLIAMS Ethmoid Forceps

Product Number	Jaw	Tip	Working Length (in)
EN0200	3 MM THRU CUT JAWS W/ BASKET, ANG UP	ANG UP	5 3/8
EN0201	ST, SERR	ST	4 3/4



WATSON-WILLIAMS Ethmoid Forceps


Product Number	Jaw Length (mm)	Tip	Size	Working Length (in)
EN0202	5.0	ST	SMALL	4 3/4
EN0203	7.0	ST	MEDIUM	4 3/4
EN0204	9.0	ST	LARGE	4 3/4



EN0202 - EN0204



EN0202



EN0203




EN0204






WILDE Nasal Forceps

Product Number	Jaw Length (mm)	Tip	Working Length (in)	Length (in)
EN0205	5 x 10	ST	4 1/2	7 1/4
EN0206	5 x 10	ANG UP	4 1/2	7 1/4



WILDE-Type Backbiter Forceps

Product Number	Jaw Size (mm)	Tip	Working Length (in)	Length (in)
EN0207	2	ST, ROTATING 360°	4 1/4	8 3/4
EN0208	2	Cd UP	3 1/8	6 7/8
EN0209	2	Cd DOWN	3 1/8	6 7/8






ASH Septum Forceps

Product Number	Tip	Length (in)
EN0210	ANG	8 3/4




COTTLE-WALSHAM Septum Forceps

Product Number	Tip	Length (in)
EN0211	ST	8 3/4



WALSHAM Septum Forceps

Product Number	Tip	Length (in)
EN0212	LEFT	8 3/4
EN0213	RIGHT	8 3/4



**Septum Morselizer Forceps,
 W/ Removable Guard**

Product Number	Length (in)
EN0214	7 7/8




EN0214



EN0214

Bone and Cartilage Crusher, W/ Clamp

Product Number
EN0215




EN0215

LANGE Fistual Hook

Product Number	Tip	Length (in)
EN0216	BLUNT	9 1/4



EN0216



EN0216

CONVERSE Skin Hooks

Product Number	Size	Length (in)
EN0217	SMALL	7 1/8
EN0218	LARGE	7 1/8



En0217 - En0218



EN0217



EN0218



EN0219



EN0219



CONVERSE Double Hook

Product Number	Blade Size (mm)	Length (in)
EN0219	7 x 12 / 8 x 12	4 3/4

MILLARD Double Hook W/ Thimble Mount

Product Number	Length (in)
EN0220	3 1/4



EN0220



EN0220

AUFRICHT Nasal Retractor

Product Number	Blade Size (mm)	Length (in)
EN0221	10 x 45	6 7/8
EN0222	10 x 45	7 3/8
EN0223	11 x 60	7 3/8



EN0221 - EN0222



EN0223

EN0223

Blair Retractor


Product Number	Tip	Blade END	Tip Size (mm)	Length (in)
EN0224	SHARP	4 PRONGS	7	5 1/4



EN0224




EN0224




CONVERSE Nasal Retractor


Product Number	Blade Size (mm)	Length (in)
EN0225	9 x 18	3 1/4
EN0226	12 x 21	3 1/4
EN0227	12 x 38	3 1/4
EN0228	12 x 53	3 1/4
EN0229	44 x 11	4 1/2



EN0225



EN0226



EN0227



EN0228




EN0229

COTTLE Nasal Retractor

Product Number	Tip	Blade ENd	Tip Size (mm)	Length (in)
EN0230	BLUNT	4 PRONGS	10	5 1/8



EN0230



EN0230

COTTLE Nasal Retractor, Double Ended

Product Number	Blade Size	Tip	Length (in)
EN0231	12 x 5 / 4 x 4.5	BALL END / BLUNT	7 7/8



EN0231



EN0231



COTTLE Retractor

Product Number	Tip	Length (in)
EN0232	BLUNT /RIGHT SHARP	5 3/4
EN0233	LEFT SHARP / BLUNT	5 3/4



EN0232 - EN0233



EN0232



EN0233

COTTLE Alar (MASING) Retractor




Product Number	Blade Size (mm)	Length (in)
EN0234	6 x 12	5 1/2
EN0235	10 x 12	5 1/2




FOMON Retractor

Product Number	Tip	Blade ENd	Tip Size (mm)	Length (in)
EN0236	BLUNT	2 PRONGS	12	5 3/4
EN0237	BALL END	2 PRONGS	13	6 1/4



GILLIES Bone and Zygoma Hook




Product Number	Tip Size (mm)	Length (in)
EN0238	9	7 1/4

PERKINS Antrum Retractor


Product Number	Right Blade	Left Blade	Length (in)
EN0239	3 PRONGS	SOLID	5
EN0240	SOLID	3 PRONGS	5




EN0239



EN0239



EN0240



EN0240



Baby WEITLANER Retractor

Product Number	Tip	Blade ENd	Length (in)
EN0241	BLUNT	3 x 3 PRONGS	4 1/4




EN0241



EN0241

WULLSTEIN Retractor

Product Number	Tip	Blade ENd	Length (in)
EN0242	BLUNT	3 x 3 PRONGS	4 1/4



JANSEN Retractor

Product Number	Tip	Blade ENd	Length (in)
EE0079	SHARP	3 x 3 PRONGS	4
EE0080	BLUNT	3 x 3 PRONGS	4
EE0081	BLUNT	3 x 4 PRONGS	4



BOIES Nasal Fracture Elevator

Product Number	Tip	Length (in)
EN0243	BLUNT	7



EN0243

GOLDMANN Septum Displacer

Product Number	Tip	Length (in)
EN0244	ST	7 7/8




EN0244

CUSHING Elevator, Narrow

Product Number	Tip Size (mm)	Tip	Length (in)
EN0245	8	SHARP	7 3/4



EN0245




EN0245

LEMPERT Elevator

Product Number	Tip Size (mm)	Blade	Length (in)
EN0246	7	Cd	7 1/4




EN0246



EN0246

LEMPERT Elevator

Product Number	Tip Size (mm)	Tip	Length (in)
EN0247	2	SLIGHT Cd	6 3/4
EN0248	3	SLIGHT Cd	6 3/4
EN0249	3	ST	7 1/2
EN0250	6	BROAD, SLIGHT Cd	7 1/4




KLEINERT-KUTZ Elevator and Dissector, Double Ended

Product Number	Tip Size (mm)	Length (in)
EN0251	2/4	7




MCKENTY Periosteal Elevator

Product Number	Tip Size (mm)	Tip	Length (in)
EN0252	4.0	Cd, SHARP	6




Ball ENd Elevator, Double Ended

Product Number	Blade Shape	Length (in)
EN0253	ANG	7 1/2




COTTLE Elevator, Double Ended

Product Number	Blade	Length (in)
EN0254	-	7 7/8
EN0255	GRADUATED	7 7/8




COTTLE Periosteal Elevator

Product Number	Tip	Width (mm)	Length (in)
EN0256	SLIGHTLY Cd	4	7
EN0257	Cd, SHARP	8	7 3/4




SEWALL Antrum Elevator

Product Number	Tip Size (mm)	Tip	Length (in)
EN0258	5.9	BLUNT	7 7/8



HAJEK-BALLENGER Septum Elevator and Dissector, Double Ended

Product Number	Tip	Length (in)
EN0259	SHARP/BLUNT	7 1/2



COTTLE (MASING) Elevator, Graduated

Product Number	Blade	Length (in)
EN0260	DOUBLE ENDED	7 7/8



PENNINGTON Elevator

Product Number	Tip	Length (in)
EN0261	SHARP/BLUNT	8 1/4



PIERCE Submucous Dissector

Product Number	Length (in)
EN0262	7 1/2




ROGERS Submucous Dissector

Product Number	Length (in)
EN0263	8 3/4




HENNER Elevator, Double Ended

Product Number	Tip Size (mm)	Length (in)
EN0264	2/2	8 1/2




HOWARTH Elevator, Double Ended

Product Number	Tip	Length (in)
EN0265	SHARP/BLUNT	8 1/2




KILLIAN Elevator, Double Ended

Product Number	Tip	Length (in)
EN0266	SHARP/BLUNT	7 1/4



Spatula and Packer, Double Ended

Product Number	Length (in)
EN0267	5 3/4




HALLE Elevator

Product Number	Tip Size (mm)	Tip	Length (in)
EN0268	3.3	BLUNT	7



EN0268





EN0268

KOENIG Elevator


Product Number	Tip Size (mm)	Tip	Length (in)
EN0269	4.5	SLIGHTLY Cd	6 1/8
EN0270	4.5	STRONGLY Cd	6 1/8



EN0269 - EN0270



EN0269



EN0270

OBWEGESER-FREER Elevator

Product Number	Tip Size (mm)	Tip	Blade Shape	Length (in)
EN0271	6	BLUNT	DOUBLE ENDED	8 1/2







EN0271



EN0271


FREER Elevator

Product Number	Tip Size (mm)	Tip	Length (in)
EN0272	4	SINGLE ENDED SHARP	7 1/2
EN0273	4	SINGLE ENDED BLUNT	7 1/2
EN0274	4 / 4	SHARP/BLUNT	7 1/4
EN0275	4.5 / 4.5	SHARP/BLUNT	7 1/2
EN0276	6.5	SHARP/BLUNT	7 1/2





FREER-YASARGIL Elevator

Product Number	Tip	Length (in)
EN0277	SHARP/SHARP	7 1/4






COTTLE Septum Speculum

Product Number	Blade Size (mm)	Blade Shape	W/ Set Screw	Length (in)
EN0278	50 x 6	TAPERED	NO	5 1/4
EN0279	76 x 7	THIN	YES	5 1/4
EN0280	90 x 7	THIN	YES	5 1/4



COTTLE Nasal Speculum

Product Number	Blade Size (mm)	Blade Shape	Size (#)	Length (in)
EN0281	20	Xtra Thin	1	5 1/4
EN0282	38	Xtra Thin	2	5 1/4
EN0283	41	Xtra Thin	3	5 1/4
EN0284	60	Xtra Thin	4	5 1/4






COTTLE Septum Speculum (continued)

Product Number	Blade Size (mm)	Blade Shape	Length (in)
EN0285	50	THIN	5 5/8



KILLIAN Septum Speculum, Aseptic

Product Number	Blade Size (mm)	Length (in)
EN0286	36 x 7	5 1/2
EN0287	56 x 7	5 1/2
EN0288	75 x 7	5 1/2
EN0289	90 x 7	5 1/2




KILLIAN Septum Speculum W/ Screw Joint

Product Number	Blade Size (mm)	Length (in)
EN0290	36 x 7	5 3/4
EN0291	56 x 7	5 3/4
EN0292	75 x 7	5 3/4
EN0293	90 x 7	5 3/4






EN0290 - EN0293



VIENNA Nasal Speculum

Product Number	Blade Size (mm)	Length (in)
EN0294	SMALL	5 1/2
EN0295	MEDIUM	5 1/2
EN0296	LARGE	5 1/2





EN0297 - EN0299

VIENNA Nasal Speculum (continued)

Product Number	Blade Size (mm)	Tip Size (mm)	Length (in)
EN0297	24	1	5 1/4
EN0298	30	3	5 1/4
EN0299	34	4	5 1/4

24 mm



EN0297

30 mm



EN0298

34 mm



EN0299

HARTMANN Nasal Speculum

Product Number	Blade Size (mm)	Length (in)
EN0300	LIGHT	5 3/4
EN0301	SMALL	6 1/4
EN0302	MEDIUM	6 1/4
EN0303	LARGE	6 1/4



EN0300

31 mm



EN0300




EN0301 - EN0303

28 mm



EN0301

31 mm



EN0302

32 mm




EN0303




EN0304 - EN0306

**HARTMANN-HALLE Nasal Speculum
W/ Screw Joint**

Product Number	Blade Size (mm)	Length (in)
EN0304	SMALL	6 1/8
EN0305	MEDIUM	6 1/8
EN0306	LARGE	6 1/8



EN0304



EN0305




EN0306

**TIECK-HALLE Nasal Speculum W/ Screw
Joint, Pediatric**


Product Number	Blade Size (mm)	Length (in)
EN0307	SMALL	5 1/4
EN0308	LARGE	5 1/4



EN0307 - EN0308




EN0307



EN0308


BECKMANN Nasal Speculum W/ Screw Joint, Pediatric

Product Number	Blade Size (mm)	Blade Shape	Length (in)
EN0309	SMALL	LAT-Cd	5 1/2
EN0310	LARGE	LAT-Cd	6



SONNENSCHN Nasal Speculum

Product Number	Blade Size (mm)	Length (in)
EN0311	MEDIUM	6



CUSHING-LANDOLT Transssphenoidal Speculum

Product Number	Blade Size (mm)	Length (in)
EN0312	70 x 15	-
EN0313	90 x 15	-
EN0314	110 x 15	-
EN0315	SPREADER FOR EN0312 - EN0314	8 1/4



PAPAVERO-HARTMANN (CASPAR) Transssphenoidal Speculum, Slim, W/ Key


Product Number	Blade Size (mm)	Description	Length (in)
EN0316	80 x 11	-	-
EN0317	90 x 13	-	-
EN0318	110 x 15	-	-
EN0319	-	HEX EY	-
EN0320	34 x 9	SPREADER FOR EN0316 - EN0318	5 1/2



EN0320

Antrum Rasp

Product Number	Tip	Length (in)
EN0321	SHARP	7
EN0322	BLUNT	7




WEINER Antrum Rasp Trocar

Product Number	Tip	Length (in)
EN0323	EXTRA SMALL, TROCAR	7 3/8




AUFRIECHT Nasal Rasp


Product Number	Tip Size (mm)	Tip	Length (in)
EN0324	9	DOWN-CUT	8 3/4
EN0325	9	UP-CUT	8 3/4



COTTLE Rasp

Product Number	Tip	Length (in)
EN0326	DOWN CUTTING	8 1/4






JOSEPH Nasal Saw, Bayonet

Product Number	Tip	Length (in)
EN0327	ST	7 3/4
EN0328	Cd RIGHT	7 3/4
EN0329	Cd LEFT	7 3/4

JOSEPH Nasal Saw

Product Number	Tip	Tip Size (mm)	Length (in)
EN0330	FINE-SERR	7.5	6 1/2
EN0331	COARSE-SERR	9	6 1/2



LEWIS Bone Rasp

Product Number	Tip	Tip Size (mm)	Length (in)
EN0332	FINE-SERR DOWN-STROKE	8	6 7/8
EN0333	COARSE-SERR DOWN-STROKE	8	6 7/8



Diamond Dust Rasp

Product Number	Tip	Length (in)
EN0334	ST	7 3/8
EN0335	Cd	7 3/8



EN0334 - EN0335



EN0334



EN0335

FOMON Bone Rasp, Double ENded

Product Number	Tip	Tip Size (mm)	Length (in)
EN0336	X-SERR FINE/COARSE	6.5 / 7.5	8 1/2
EN0337	FINE/COARSE	6.5 / 7.5	8 1/2
EN0338	FLAT, CONVEX FINE, COARSE	8/8	8 1/2
EN0339	COARSE / XTRA-COARSE	7.5 / 7.5	8 3/8



EN0336



EN0336



EN0337



EN0337



EN0338



EN0338



EN0339



EN0339

GALLAHER File Rasp, Coarse

Product Number	Tip	Length (in)
EN0340	6.5	7 1/2



EN0340



EN0340



EN0341



EN0341

MALTZ-ANDERSON Nasal Rasp

Product Number	Tip Size (mm)	Length (in)
EN0341	5.7	6 1/2

MALTZ Bone Rasp

Product Number	Tip	Tip Size (mm)	Length (in)
EN0342	DOWN STROKE	9	6 1/2
EN0343	DOWN STROKE	33	6 1/2
EN0344	UP STROKE	33	6 1/2



EN0342



EN0342



EN0343



EN0343



EN0344



EN0344



EN0345



EN0345

PARKES Nasal Rasp

Product Number	Tip Size (mm)	Length (in)
EN0345	7	7 7/8

BALLENGER Chisel, V-shaped

Product Number	Tip Size (mm)	Length (in)
EN0346	2	5 7/8
EN0347	4	5 7/8
EN0348	8	5 7/8



BALLENGER Septum Gouge

Product Number	Tip Size (mm)	Length (in)
EN0349	6	6


CONVERSE Guarded Chisel

Product Number	Tip Size (mm)	Blade Shape	Length (in)
EN0350	6	Cd RIGHT	7 1/2
EN0351	6	Cd LEFT	7 1/2




COTTLE Chisel W/ Depth Marks

Product Number	Tip Size (mm)	Length (in)
EN0352	4	7
EN0353	7	7
EN0354	9	7
EN0355	12	7



FREER Chisel



Product Number	Tip Size (mm)	Blade Shape	Length (in)
EN0356	4	ST	6 1/4
EN0357	4	Cd	6 1/4

PAPAVERO Bayonet Chisel

Product Number	Tip Size (mm)	Length (in)
EN0358	4	8 3/4



FREER Septum Gouge



Product Number	Tip Size (mm)	Length (in)
EN0359	4	6


KILLIAN Septum Gouge, Bayonet



Product Number	Tip Size (mm)	Length (in)
EN0360	6	6 3/4

PARKES Hump Gouge, Double Guard

Product Number	Tip Size (mm)	Blade Shape	Length (in)
EN0361	6	ST	7
EN0362	10	ST	7



ANDERSON-NEIVERT Osteotome

Product Number	Tip Size (mm)	Blade Shape	Length (in)
EN0363	5	Cd RIGHT	8
EN0364	5	Cd LEFT	8



EN0363




EN0364

CHERRY Osteotome

Product Number	Tip Size (mm)	Length (in)
EN0365	3	8
EN0366	4	8
EN0367	5	8
EN0368	6	8
EN0369	7	8



EN0365 - EN0369



EN0365 - EN0369

CINELLI Osteotome, Double Guard

Product Number	Tip Size (mm)	Blade Shape	Length (in)
EN0370	10	Cd	6 1/8



EN0370

CINELLI Osteotome, Double Guard

Product Number	Tip Size (mm)	Blade	Length (in)
EN0371	10	ST	6 1/4
EN0372	12	ST	6 1/4
EN0373	14	ST	6 1/4
EN0374	16	ST	6 1/4



EN0371 - EN0374



EN0371



EN0372



EN0373



EN0374



EN0375 - EN0380



EN0375 - EN0380

CONVERSE Osteotome


Product Number	Tip Size (mm)	Length (in)
EN0375	2	7
EN0376	4	7
EN0377	6	7
EN0378	8	7
EN0379	10	7
EN0380	12	7

COTTLE Osteotome, W/ Depth Mark

Product Number	Tip Size (mm)	Length (in)
EN0381	4	7
EN0382	6	7
EN0383	8	7
EN0384	10	7
EN0385	12	7




EN0380 - EN0385



EN0380 - EN0385

COTTLE Osteotome, Curved

Product Number	Tip Size (mm)	Length (in)
EN0386	6	7




COTTLE Osteotome W/ Cross Bar



Product Number	Tip Size (mm)	Blade	Length (in)
EN0387	2	ST	7
EN0388	4	ST	7
EN0389	6	ST	7
EN0390	7	ST	7
EN0391	9	ST	7
EN0392	12	ST	7

TESSIER Osteotome

Product Number	Tip Size (mm)	Blade	Length (in)
EN0393	2.5	ST	6 1/4
EN0394	3.5	ST	6 1/4
EN0395	5	ST	6 1/4
EN0396	7	ST	6 1/4
EN0397	10	ST	6 1/4
EN0398	14	ST	6 1/4
EN0399	18	ST	6 1/4
EN0400	5	Cd	6 1/4
EN0401	7	Cd	6 1/4
EN0402	10	Cd	6 1/4



DAUTREY Osteotome

Product Number	Tip Size (mm)	Blade	Length (in)
EN0403	4	ST	6 3/4
EN0404	6	ST	6 3/4
EN0405	6	SLIGHTLY Cd	6 3/4
EN0406	6	STRONGLY Cd	6 3/4
EN0407	8	ST	6 3/4
EN0408	8	SLIGHTLY Cd	6 3/4
EN0409	8	STRONGLY Cd	6 3/4




EN0403




EN0403




EN0404 - EN0406



EN0404



EN0405



EN0406



EN0407 - EN0409



EN0407



EN0408




EN0409


FOMON Osteotome, Single Guard



EN0410 - EN0411



EN0410



EN0411

Product Number	Tip Size (mm)	Length (in)
EN0410	5	6 3/4
EN0411	8	6 3/4

MASING Osteotome, Single Guard

Product Number	Blade	Length (in)
EN0412	ST	6 7/8
EN0413	LEFT Cd	6 7/8
EN0414	RIGHT Cd	6 7/8



EN0412



EN0413



EN0414

OBEWESGESER Nasal Osteotome



EN0415



EN0416



EN0417

Product Number	Tip Size (mm)	Length (in)
EN0415	6	7 1/4
EN0416	8	7 1/4
EN0417	4	7 1/2



Osteotome W/ Hexagonal Handle



Product Number	Tip Size (mm)	Length (in)
EN0418	2	6 1/4
EN0419	3	6 1/4
EN0420	4	6 1/4
EN0421	5	6 1/4
EN0422	6	6 1/4
EN0423	7	6 1/4
EN0424	8	6 1/4
EN0425	9	6 1/4
EN0426	10	6 1/4
EN0427	11	6 1/4
EN0428	12	6 1/4


RUBIN Osteotome, W/ Finger Grip

Product Number	Tip Size (mm)	Length (in)
EN0429	14	6 1/4
EN0430	10	6 1/2
EN0431	14	6 1/2



Frontal Ostium Seeker Probe, Double Ended

Product Number	Tip	Length (in)
EN0432	Cd RIGHT, ANG	8 1/4



Maxillary Sinus Seeker Probe, Double Ended


Product Number	Tip	Tip Size (mm)	Length (in)
EN0433	BALL	1.5/ 2	7 1/2



SKILLERN Sinus Probe

Product Number	Tip	Length (in)
EN0434	MALLEABLE W/ BEADS AT 50, 70 AND 90 MM	6 3/4






BALLENGER Swivel Knife, Bayonet

Product Number	Blade Size (mm)	Length (in)
EN0435	4 x 3	8 1/4
EN0436	5 x 4	8 1/4
EN0437	6 x 5	8 1/4

BALLENGER Swivel Knife, Straight

Product Number	Blade Width (mm)	Length (in)
EN0438	3	7 1/2
EN0439	3	7 3/8
EN0440	4	7 3/8
EN0441	5	7 1/2




CONVERSE Septum Knife, Button END

Product Number	Blade	Length (in)
EN0442	Cd	6

COTTLE Septum Knife, Round ENd

Product Number	Blade	Length (in)
EN0443	ST	5 1/2



FREER Septum Knife

Product Number	Blade Shape	Blade Size (mm)	Length (in)
EN0444	"D" SHAPED	3	6
EN0445	"D" SHAPED	3	6 1/8
EN0446	ANG	5 x 8	7



EN0444



EN0444



EN0445



EN0445



EN0446

JOSEPH Knife, Button End

Product Number	Tip	Length (in)
EN0447	ST	5 1/2
EN0448	Cd	6 7/8



EN0447



EN0447



EN0448



EN0448

JOSEPH Knife, Double Edge

Product Number	Blade Shape	Blade Size (mm)	Length (in)
EN0449	ST	5 x 23	5 3/4
EN0450	Cd	5 x 23	5 3/4



EN0449 - EN0450



EN0449



EN0450



EN0450

Suction Tube, Angled

Product Number	Tip	Tip Size (mm)	Length (in)
EN0451	ANG	13	3



Nasal Suction Tube W/ Luer Hub, Tapered

Product Number	Tip Size (mm)	Length (in)
GS1549	5	7




FERGUSON-FRAZIER Suction Tube

Product Number	Tip	Tip Size (FR)	Tip Size (mm)	Length (in)
EN0452	ANG	8	7 7/8	10 1/4
EN0453	ANG	10	7 7/8	10 1/4
EN0454	ANG	14	4	6 1/2
EN0455	ANG	14	5 7/8	8



FRAZIER Suction Tube

Product Number	Tip Size (FR)	Description
GS1605	7	-
GS1606	7	SUCTION STYLET W/OLIVE TIP
GS1607	8	SUCTION STYLET W/OLIVE TIP
GS1608	9	-
GS1609	-	SUCTION STYLET W/OLIVE TIP
GS1610	10	-
GS1611	10	SUCTION STYLET W/OLIVE TIP
GS1612	11	-
GS1613	11	SUCTION STYLET W/OLIVE TIP
GS1614	12	-
GS1615	12	SUCTION STYLET W/OLIVE TIP




FRAZIER Suction Tube, Straight

Product Number	Tip	Working Length (in)	Length (in)
EN0456	ST	5 1/2	6 7/8



FRAZIER Suction Tube W/Luer Hub

Product Number	Tip Size (FR)	Length (in)
GS1628	9	7
GS1629	15	7



EICKEN Antrum Suction Tube, Bulbous

Product Number	Tip	Diameter (mm)	Length (in)
EN0457	SMALL	2.5	5 1/8
EN0458	LARGE	3	5 1/8
EN0459	LONG CURVE	3	5 5/8
EN0460	SHORT CURVE	3	5 5/8
EN0461	SHORT CURVE	4	5 5/8

KILLIAN Antrum Cannula

Product Number	Tip	Tip Size (#)	Length (in)
EN0462	SMALL	1	5
EN0463	LARGE	2	5




SACHS Suction Tube

Product Number	Tip Size	Length (in)
GS1707	SMALL	8
GS1708	LARGE	8



GORNEY Suction Dissector, Elevator

Product Number	Tip Size (mm)	Length (in)
EN0464	6-10	7
EN0465	3	8 1/2



HAJEK Irrigating Cannula

Product Number	Diameter (mm)	Length (in)
EN0466	3.5	6 1/2

HOUSE Cutoff Adapter

Product Number	Tip	Tip Size (mm)
EN0467	TUBING END	4



Stopcock Luer Adapter, Female to Male

Product Number	Description
EN0468	FEMALE TO MALE LUER LOCK W/ FEMALE RIGHT

Luer Adapter

Product Number	Description
EN0469	MALE LUER LOCK- TO - HOSE





COTTLE Mallet, Flat/Round, 235gr

Product Number	Diameter (mm)	Length (in)
EN0470	30	7 1/4



EN0470

CRANE Mallet

Product Number	Diameter (mm)	Length (in)
EN0471	32	6 1/4



EN0471

Mallet, 155gr

Product Number	Diameter (mm)	Length (in)
EN0472	25	7 1/2



EN0472

Lead Mallet

Product Number	Diameter (mm)	Handle	Description	Length (in)
EN0473	25	LARGE	LEAD FILLED	7 1/16



EN0473

Mallet W/ Removable Nylon Discs, 135gr

Product Number	Diameter (mm)	Description	Length (in)
EN0474	25	COMPLETE	7 1/4
EN0475	25	NYLON SPARE DISCS, 2 SET	-



EN0474 - EN0475

QUISLING Intranasal Hammer

Product Number	Length (in)
EN0476	6 3/4



EN0476



EN0477 - EN0480

○
EN0477

○
EN0478

○
EN0479

○
EN0480

Bone Tamp

Product Number	Diameter (mm)	Length (in)
EN0477	3	6 1/4
EN0478	5	6 1/4
EN0479	7	6 1/4
EN0480	9	6 1/4

Antrum Curette

Product Number	Tip	Tip Size (mm)	Length (in)
EN0481	OBLONG	-	7 1/2
EN0482	OVAL CUP	7 x 11	7 1/2
EN0483	OVAL CUP	8 x 11	7 1/2
EN0484	OVAL CUP	9 x 13	7 1/2
EN0485	OVAL CUP, ANG	9.5 x 4	7 1/2
EN0486	OVAL CUP, ANG	6 x 3	7 1/2
EN0487	55° Cd	6 x 3	7 1/2
EN0488	90° Cd	6 X3	7 1/2
EN0489	BACKWARD CUTTING	7 X10	7 1/2
EN0490	FORWARD CUTTING	7 X10	7 1/2



EN0481 - EN0484



EN0481



EN0483



EN0482




EN0484



EN0485 - EN0488



EN0485



EN0487



EN0486



EN0488



EN0489



EN0489



EN0490

EN0490

BARTH Bone Curette, Double Ended



Product Number	Tip Size (mm)	Length (in)
EN0491	2.3 x 6.0 / 2.7 x 7.0	7

COAKLEY Antrum Curette

Product Number	Blade	Size (#)	Length (in)
EN0492	-	SET OF 6 INCLUDES EN0493 - EN0498	-
EN0493	SLIGHT ANG	1	6 3/4
EN0494	45°	2	6 3/4
EN0495	REV ANG	3	6 3/4
EN0496	SLIGHT ANG	4	6 3/4
EN0497	45°	5	6 3/4
EN0498	45°	6	6 3/4



EN0493 - EN0498



EN0493



EN0496



EN0494



EN0497



EN0495



EN0498



EN0499



EN0499

JACOBSON Bone Curette, Double Ended

Product Number	Tip Size (mm)	Length (in)
EN0499	4.3 x 6.5 / 4.5 x 8.5	5 1/2

JANSEN Bone Curette, Double Ended

Product Number	Tip Size (mm)	Length (in)
EN0500	4.3 / 3.6	6 1/4




EN0500



EN0500



EN0501



EN0501



EN0502



EN0502

VOLKMANN Bone Curette, Double Ended

Product Number	Tip Size (mm)	Length (in)
EN0501	7.5 x 21.0 / 6.0 x 11.0	5 1/8
EN0502	7.5 x 21.0 / 6.0 x 11.0	6 3/4

WILLIGER Bone Curette, Double Ended

Product Number	Tip Size (mm)	Length (in)
EN0503	2.0 x 3.5 / 3.0 x 4.5	5 3/4




EN0503



EN0503



EN0504 - EN0507





LEMPERT Bone Curette

Product Number	Size (#)	Length (in)
EN0504	000000	8
EN0505	00000	8
EN0506	0000	8
EN0507	2	8

SPRATT Mastoid Curettes

Product Number	Tip Size (mm)	Size (#)	Length (in)
EN0508	2.5 x 4	00	6 1/2
EN0509	3 x 5	0	6 1/2
EN0510	3.5 x 6	1	6 1/2
EN0511	4.5 x 8	3	6 1/2
EN0512	5.5 x 9	4	6 1/2







BLUMENTHAL Bone Rongeur

Product Number	Tip Size (mm)	Size (#)	Length (in)
EN0513	3.5 x 25	Cd	6 1/8

BOEHLER Bone Rongeur, Double Action

Product Number	Jaw Size (mm)	Jaw	Jaw ENd	Length (in)
EN0514	3.0 x 27	Cd	NARR	6 1/8
EN0515	3.4 x 30	Cd	WIDE	6 1/8




Bone Rongeur

Product Number	Jaw Size (mm)	Jaw	Length (in)
EN0516	3.8 x 19	Cd	5 1/2

FRIEDMANN Micro Rongeur

Product Number	Jaw Size (mm)	Jaw	Length (in)
EN0517	1.3 x 19	ST	5 1/2
EN0518	1.3 x 20	Cd	5 1/2



JANSEN Bone Rongeur

Product Number	Jaw Size (mm)	Jaw	Length (in)
EN0519	4.2 x 35	ST	6 7/8
EN0520	4.5 x 20	Cd	6 7/8
EN0521	3.8 x 37	ANG	7 1/2



EN0519 - EN0520



EN0519



EN0520




EN0521



EN0521

MEAD Bone Rongeur

Product Number	Jaw Size (mm)	Jaw	Length (in)
EN0522	3.5 x 30	Cd	6 1/2



ZAUFAL-JANSEN Bone Rongeur,
Double Action

Product Number	Jaw Size (mm)	Jaw	Length (in)
EN0523	4.5 x 29	Cd	7




ALEXANDER Bone Rongeur

Product Number	Jaw Size (mm)	Jaw	Length (in)
EN0524	3.2 x 27	ANG	5 1/2



BANE-HARTMANN Bone Rongeur

Product Number	Jaw Size (mm)	Jaw	Length (in)
EN0525	4.0 x 32	Cd	6 7/8




BEYER Bone Rongeur

Product Number	Jaw Size (mm)	Jaw	Handle	Length (in)
EN0526	3.2 x 32	Cd	Single Action	7
EN0527	4.0 x 30	Cd	Double Action	7



LEMPERT Bone Rongeur




Product Number	Jaw Size (mm)	Jaw	Length (in)
EN0528	2.0 x 23	DEL, ST	6 1/2
EN0529	2.2 x 27	ST	7 7/8
EN0530	2.5 x 30	Cd	6 1/4
EN0531	2.2 x 27	Cd	7 7/8



LUER Bone Rongeur

Product Number	Jaw Size (mm)	Jaw	Length (in)
EN0532	4.2 x 20	Cd	6
EN0533	4.2 x 26	STRONG Cd	6 1/8



LUER-FRIEDMANN Bone Rongeur

Product Number	Jaw Size (mm)	Jaw	Length (in)
EN0534	2.2 x 22	Cd	5 3/4



Scan to know more

Thomson Surgicals Pvt. Ltd

31/26, 1st Floor, 6th Cross, Bhavani Nagar,
Bengaluru - 560 029, Karnataka, India

P +91 80 26685076 M +91 94808 12041

E mail@thomsonsurgicals.com

Mfg Units Bengaluru - Karnataka | Palakad - Kerala

W www.thomsonsurgicals.com

