



Obstetrics /Gynecology

“Empowering Precision in Every Incision.”





Thomson surgicals instruments possesses several key features that contribute to its effectiveness in various medical procedures:

Precision and Accuracy:

Instruments designed for precise and accurate performance to enhance surgical outcomes.

Ergonomic Design:

Ergonomically designed handles for comfortable grip and reduced hand fatigue during long surgeries.

Durability and Material Quality:

Use of high-quality materials such as stainless steel grade 420, 410, 316, and 304 or titanium medical grade, Copper, Brass and Hylam Material for durability and resistance to corrosion.

Customization:

Offering customizable options based on surgeon preferences or specific surgical procedures.

Innovative Coatings:

Integration of advanced coatings for reduced friction, enhanced visibility, long life and ease of cleaning.

Minimally Invasive Solutions:

Developing instruments suitable for minimally invasive procedures, promoting smaller incisions and faster recovery times.


Modularity:

Modular instrument systems that allow for flexibility and adaptability in various surgical settings.

Easy Maintenance and Sterilization:

Designing instruments that are easy to clean, sterilize, and maintain to ensure optimal performance and infection control.

The combination of these features ensures that the Thomson surgicals instruments contribute positively to the surgical process, providing optimal exposure, accessibility, and safety during procedures.



Content

1. Scissors	01 - 09
2. Needle Holders	10
3. Forceps	11 - 48
4. Clamps	49 - 58
5. Retractors	59 - 60
6. Vaginal Retractors and Specula	61 - 79
7. Endospecula	80
8. Dilators and Probes	81 - 86
9. Curettes	87 - 95
10. Accessories	96
11. Numerical Index	97 - 102

SIMS Scissors

Product Number	Tip	Blade	Length (in)
OG0001	S/S	ST	8
OG0002	S/S	Cd	8



OG0001




OG0002



OG0001

SIMS Scissors

Product Number	Tip	Blade	Length (in)
OG0003	B/B	ST	8
OG0004	B/B	ST	9
OG0005	B/B	Cd	8
OG0006	B/B	Cd	9



OG0003 - OG0004



OG0003



OG0004



OG0005



OG0006




SIMS Scissors (continued)

Product Number	Tip	Blade	Length (in)
OG0007	S/B	ST	8
OG0008	S/B	Cd	8
OG0009	S/B	Cd	9



OG0007



OG0007



OG0008 - OG0009

WERTHEIM Scissors

Product Number	Tip	Blade	Length (in)
OG0010	B/B	Cd	6
OG0011	B/B	Cd	8
OG0012	B/B	Cd	9



OG00010



OG00011



OG00012



OG0010 - OG0012



OG0013

JORGENSEN Scissors


Product Number	Tip	Blade	Length (in)
OG0013	B/B	Cd	9



OG0013

SIEBOLD Uterine Scissors

Product Number	Tip	Blade	Length (in)
OG0014	B/B	Cd	10



OG0014




OG0014



OG0015

BOZEMANN Uterine Scissors

Product Number	Tip	Blade	Length (in)
OG0015	B/B	Cd	9



OG0015

SIMS-SIEBOLD Uterine Scissors

Product Number	Tip	Blade	Length (in)
OG0016	B/B	ST	9



OG0016



OG0016



OG0017 - OG0018

Umbilical Cord Scissors

Product Number	Tip	Blade	Length (in)
OG0017	B/B	LAT-Cd	5
OG0018	B/B	LAT-Cd	6



OG0017



OG0018



OG0017

Umbilical Cord Scissors

Product Number	Length (in)
OG0019	4



OG0019



OG0019



OG0020

SCHUHMACHER Umbilical Cord Scissors

Product Number	Length (in)
OG0020	6




OG0020

BRAUN-STADLER Episiotomy Scissors

Product Number	Blade	Length (in)
OG0021	SERR	6
OG0022	SERR	8



OG0021 - OG0022



OG0021



OG0022



OG0023 - OG0024

SIMS Uterine Wavecut Scissors with tungsten carbide edges

Product Number	Tip	Blade	Length (in)
OG0023	B/B	Cd	8
OG0024	B/B	Cd	9



OG0023



OG0024

METZENBAUM Wavecut Dissecting Scissors, Heavy with tungsten carbide edges

Product Number	Tip	Blade	Length (in)
GS0214	B/B	Cd	9



GS0214




GS0214



MAYO Dissecting Scissors with tungsten carbide edges

Product Number	Tip	Blade	Length (in)
GS0176	B/B	ST	9
GS0177	B/B	Cd	9
GS0178	B/B	ST	14
GS0179	B/B	Cd	14




GS0176




GS0176



GS0177



GS0178



GS0179

SIMS Uterine Scissors with tungsten carbide edges

Product Number	Tip	Blade	Length (in)
OG0025	B/B	Cd	8
OG0026	B/B	Cd	9




OG0025



OG0025 - OG0026

WERTHEIM Scissors with tungsten carbide edges

Product Number	Length (in)
OG0027	9



OG0027



OG0028 - OG0034

Uterine Scissors with tungsten carbide edges


Product Number	Tip	Length (in)
OG0028	SLIGHT Cd	9
OG0029	STRONG Cd	9
OG0030	ANG	9
OG0031	ANG	10
OG0032	SLIGHT Cd	12
OG0033	STRONG Cd	12
OG0034	ANG	12



OG0028




OG0032




OG0029



OG0033



OG0030



OG0034



OG0031



GS0282

HEANEY Needle Holder

Product Number	Tip	Blade	Length (in)
GS0282	Cd	SERR	8




GS0282

HEANEY Needle Holder with Tungsten Carbide Inserts

Product Number	Tip	Blade	Length (in)
GS0350	Cd	SERR	8
GS0351	Cd	SERR	10



0.5
GS0350



0.5
GS0351




GS0350 - GS0351

MAIER Forceps

Product Number	Jaw	Length (in)
OG0035	ST	10
OG0036	Cd	10



OG0035



OG0036




OG0035 - OG0036

MAIER Forceps, W/ Ratchet

Product Number	Jaw	Length (in)
OG0037	ST	10
OG0038	Cd	10



OG0037



OG0038



OG0037 OG0038




SIMS-MAIER Sponge Forceps, W/ Ratchet

Product Number	Jaw	Length (in)
OG0039	ST	11
OG0040	Cd	11



OG0039



OG0040



OG0039 - OG0040

PELKMAN Forceps

Product Number	Tip	Jaw	Length (in)
OG0041	SERR	ST	10
OG0042	SERR	Cd	10



OG0041 - OG0042



OG0041




OG0042

BOZEMANN Forceps


Product Number	Shaft	Jaw	Length (in)
OG0043	ST	SERR	10
OG0044	S-Cd	SMOOTH	10
OG0045	S-Cd	SERR	10
OG0046	S-Cd	W/ GROOVE	10



OG0043



OG0043



OG0045




OG0044 - OG0046



BOZEMANN-DOUGLAS Forceps


Product Number	Tip	Length (in)
OG0047	SERR	10



OG0047

ADAIR Uterine Forceps

Product Number	Teeth	Length (in)
OG0048	1 x 1	10




OG0048



OG0048

POZZI Tenaculum Forceps, Standard

Product Number	Teeth	Length (in)
OG0049	1 x 1	7



OG0049



OG0049



OG0050

BARRETT Uterine Tenaculum Forceps

Product Number	Teeth	Length (in)
OG0050	1 x 1	7




OG0050



BRAUN Tenaculum Forceps

Product Number	Teeth	Length (in)
OG0051	1 x 1	10



OG0051



OG0051

SCHROEDER Tenaculum Forceps

Product Number	Shaft	Teeth	Length (in)
OG0052	ST	1 x 1	10



OG0052



OG0052

SCHROEDER Tenaculum Forceps

Product Number	Shaft	Teeth	Length (in)
OG0053	ST	2 x 2	10
OG0054	Cd	2 x 2	10
OG0055	ST	3 x 3	11



OG0053



OG0054



OG0055



OG0055



OG0053



OG0056



OG0056

POZZI Tenaculum Forceps, Heavy


Product Number	Teeth	Length (in)
OG0056	1 x 1	10

DUPLAY Tenaculum Forceps, Delicate

Product Number	Teeth	Length (in)
OG0057	1 x 1	10



OG0057



OG0057

ABEL Tenaculum Forceps

Product Number	Teeth	Length (in)
OG0058	1 x 1	11



OG0058



OG0058



OG0059

HENROTIN Uterine Vulsellum Forceps

Product Number	Teeth	Length (in)
OG0059	2 x 3	8



OG0059





OG0060

GORDON Uterine Vulsellum Forceps

Product Number	Teeth	Length (in)
OG0060	3 x 3	8

CZERNY Uterine Tenaculum Forceps

Product Number	Teeth	Length (in)
OG0061	4 x 4	8



OG0061

OG0061



OG0062

OG0062

PRATT Tenaculum Forceps

Product Number	Teeth	Length (in)
OG0062	4 x 4	10

PRATT Tenaculum Forceps


Product Number	Teeth	Length (in)
OG0063	4 x 4	11




OG0063



OG0063



OG0064



OG0064

TEALE Tenaculum Forceps

Product Number	Teeth	Length (in)
OG0064	3 x 4	10



JACOBS Vulsellum Forceps


Product Number	Teeth	Length (in)
OG0065	2 x 2	8



OG0065

LITTLEWOOD Tenaculum Forceps

Product Number	Teeth	Length (in)
OG0066	2 x 3	7



OG0066

MUSEUX Tenaculum Forceps, Small


Product Number	Teeth	Length (in)
OG0067	2 x 2	8




OG0067

MUSEUX Tenaculum Forceps


Product Number	Jaw	Length (in)
OG0068	ST	9
OG0069	ST	9
OG0070	ST	9



OG0068



OG0069





OG0070



OG0068 - OG0070

MUSEUX Tenaculum Forceps

Product Number	Jaw	Length (in)
OG0071	Cd	9




OG0071




OG0072 - OG0073

MUSEUX Tenaculum Forceps

Product Number	Jaw	Size	Length (in)
OG0072	ST	LARGE	10
OG0073	ST	X-LARGE	10



OG0072



OG0073

FLETCHER-VAN DOREN Uterine Polyp, Sponge Holding Forceps

Product Number	Length (in)
OG0074	9



OG0074

FLETCHER Sponge Forceps

Product Number	Tip	Length (in)
OG0075	ST	10



OG0075



OG0075

OVERSTREET Endometrial Polyp Forceps

Product Number	Jaw	Length (in)
OG0076	HOLLOW, FENESTRATED	10



OG0076



OG0076



HEYWOOD-SMITH Polyp Forceps

Product Number	Length (in)
OG0077	9



OG0077



OG0077



OG0078



OG0078

NOTO Polyp Forceps

Product Number	Length (in)
OG0078	9

DOYEN Ovum Forceps

Product Number	Length (in)
OG0079	9




OG0079



OG0079



OG0080



OG0080

Atraumatic Ovum Forceps

Product Number	Tip	Length (in)
OG0080	SMALL	10



COLLIN Uterine Elevating Forceps

Product Number	Length (in)
OG0081	11



OG0081




OG0082

COLLIN Uterine Elevating Forceps, Movable Jaw

Product Number	Blade	Length (in)
OG0082	ST	10

DARTIGUES Elevating Forceps

Product Number	Length (in)
OG0083	10



SOMER Uterine Elevating Forceps

Product Number	Length (in)
OG0084	9



HULKA-KENNWICK Uterine Elevating Forceps

Product Number	Length (in)
OG0085	11



OG0085



OG0085

IOWA Membrane Puncturing Forceps

Product Number	Length (in)
OG0086	10





OG0086



OG0086

DOYEN Myoma Screw


Product Number	Length (in)
OG0087	6
OG0088	7





NAEGELE Obstetrical Forceps

Product Number	Length (in)
OG0089	14



OG0089



OG0090

KIELLAND Obstetrical Forceps

Product Number	Length (in)
OG0090	16

Obstetrical Forceps

Product Number	Length (in)
OG0091	12
OG0092	14



OG0091



OG0092

DE LEE Obstetrical Forceps

Product Number	Length (in)
OG0093	12
OG0094	14



OG0093 - OG0094



OG0095

ELLIOTT Obstetrical Forceps

Product Number	Length (in)
OG0095	15

PIPER Obstetrical Forceps, F/ Breach

Product Number	Length (in)
OG0096	17



OG0096



LUIKART-SIMPSON Obstetrical Forceps



Product Number	Length (in)
OG0097	14

OG0097

LUIKART-MCLANE Obstetrical Forceps

Product Number	Length (in)
OG0098	14



OG0098

LUIKART-KIELLAND Obstetrical Forceps



Product Number	Length (in)
OG0099	15

OG0099

SIMPSON C-Section Forceps

Product Number	Length (in)
OG0100	9



OG0100

BARNES-NEVILLE Obstetrical Forceps, With Traction Handle

Product Number	Length (in)
OG0101	17



OG0101



OG0102

BILL Traction Handle

Product Number	Description
OG0102	FOR W/STANDARD OB FORCEPS



LUIKHART Obstetrical Forceps

Product Number	Length (in)
OG0103	16



MCLANE-TUCKER Obstetrical Forceps

Product Number	Length (in)
OG0104	16



MURLESS Head Extractor

Product Number	Length (in)
OG0105	10





OG0106

SELLHEIM Obstetric Lever

Product Number	Length (in)
OG0106	16

**GREEN-ARMYTAGE Forceps,
F/ C-Section**

Product Number	Length (in)
OG0107	9



OG0107



OG0107



BARRETT Placenta Forceps

Product Number	Blade	Length (in)
OG0108	FENESTRATED	10



OG0108



OG0108

BIERER Ovum Forceps

Product Number	Jaw(mm)	Ratchet	Length (in)
OG0109	16	NO	11
OG0110	16	YES	11
OG0111	19	NO	11
OG0112	19	YES	11



OG0109, OG0111



OG0110, OG0112



WINTER Placenta and Ovum Forceps,W/O Ratchet

Product Number	Size (#)	Jaw	Length (in)
OG0113	1	ST	11
OG0114	2	ST	11



OG0113



OG0114



OG0113 - OG0114




WINTER Placenta and Ovum Forceps, W/O Ratchet (continued)

Product Number	Size (#)	Jaw	Length (in)
OG0115	1	Cd	11
OG0116	2	Cd	11
OG0117	3	Cd	11



OG0115



OG0116



OG0117



OG0115 -OG0117

WINTER Placenta and Ovum Forceps, W/ Ratchet

Product Number	Size(#)	Jaw	Length (in)
OG0118	2	Cd	11
OG0119	3	Cd	11



OG0118



OG0119



OG0118 - OG0119



**SAENGER Placenta and Ovum Forceps,
W/ Ratchet**

Product Number	Jaw	Length (in)
OG0120	ST	11
OG0121	Cd	10
OG0122	Cd	11



OG0121



OG0120



OG0121



OG0120



OG0122



OG0122

**MCCLINTOCK Placenta and Ovum
Forceps**

Product Number	Length (in)
OG0123	9



OG0123



OG0123



**THOMS-GAYLOR Biopsy Forceps,
S-Shape**

Product Number	Length (in)
OG0124	9



OG0124



OG0124



OG0125 - OG0126

SCHUMACHER Biopsy Forceps, S-Shape

Product Number	Tip	Length (in)
OG0125	SMALL	9
OG0126	LARGE	9



OG0125

OG0126

SCHUMACHER Biopsy Rongeur

Product Number	Tip	Length (in)
OG0127	SMALL	9
OG0128	LARGE	9



OG0127

OG0128

OG0127 - OG0128




BERGER Biopsy Forceps

Product Number	Length (in)
OG0129	8




TISCHLER Cervical Biopsy Punch Forceps

Product Number	Jaw Size(mm)	Length (in)
OG0130	5 x 6 x 3	8
OG0131	7 x 3	11



TISCHLER Biopsy Alligator Forceps

Product Number	Jaw Size(mm)	Jaw	Length (in)
OG0132	7 x 3	OBLONG CUP SINGLE TOOTH	10



KEVORKIAN-YOUNGE Biopsy Forceps


Product Number	Jaw Size (mm)	Working Length (in)
OG0133	8.5 x 3.5	10





WITTNER Biopsy Forceps

Product Number	Jaw Size(mm)	Jaw	Working Length (in)
OG0134	4 x 8	ST, JAWS W/TEETH	9



OG0134

OG0134

EPPENDORF Biopsy Forceps

Product Number	Jaw Size(mm)	Working Length (in)
OG0135	8 x 4	8



OG0135

OG0135



OG0136 - OG0137

PHANEUF Hysterectomy Forceps

Product Number	Jaw	Length (in)
OG0136	ST	8
OG0137	ANG	8



OG0136



OG0137

HOLZBACH Hysterectomy Forceps

Product Number	Jaw	Length (in)
OG0138	ST	10
OG0139	Cd	10



OG0138




OG0139



OG0138 - OG0139

HEANEY-REZEK Hysterectomy Forceps

Product Number	Jaw	Length (in)
OG0140	ST	8



OG0140



OG0140



HEANEY-BALLENTINE Hysterectomy Forceps, Longitudinal Grooves

Product Number	Jaw	Length (in)
OG0141	ST	8
OG0142	ST	10
OG0143	ST	12
OG0144	ST	14
OG0145	Cd	8
OG0146	Cd	10
OG0147	Cd	12
OG0148	Cd	14



OG0141 - OG0148



OG0141 - OG0144



OG0145 - OG0148

HEANEY Hysterectomy Forceps

Product Number	Jaw	Teeth	Length (in)
OG0149	Cd	1	8
OG0150	Cd	2	8
OG0151	Cd	2	10



OG0149



OG0149



OG0150



OG0151



OG0150 - OG0151

HEANEY Hysterectomy Forceps (continued)

Product Number	Jaw	Teeth	Length (in)
OG0152	ST	2	12
OG0153	Cd	2	12



OG0152




OG0153



OG0152 - OG0153



OG0154



OG0154

HEANEY Hysterectomy Forceps, Heavy

Product Number	Jaw	Length (in)
OG0154	Cd	8



HEANEY Hysterectomy Forceps, W/O Teeth

Product Number	Jaw	Length (in)
OG0155	ST	9
OG0156	STRONG - Cd	8
OG0157	SLIGHT - Cd	9



OG0155 - OG0157



OG0155



OG0156



OG0157

WIENER Hysterectomy Forceps, W/ Grooves

Product Number	Jaw	Tip	Length (in)
OG0158	Cd	SMOOTH	9



OG0158



OG0158

WERTHEIM Atraumatic Parametrium Clamp

Product Number	Jaw	Length (in)
OG0159	ST	10
OG0160	SLIGHT-Cd	10
OG0161	Cd	9
OG0162	UP-Cd	9



OG0159



OG0160



OG0161



OG0162




OG0159 - OG0162




WERTHEIM Parametrium Clamp

Product Number	Jaw	Length (in)
OG0163	Cd	10



OG0163



OG0163

WERTHEIM Hysterectomy Forceps

Product Number	Jaw	Length (in)
OG0164	Cd	10




OG0164



OG0164

WERTHEIM-CULLEN Clamp

Product Number	Jaw	Length (in)
OG0165	RIGHT-ANG	8



WERTHEIM Atraumatic Clamp

Product Number	Jaw	Length (in)
OG0166	RIGHT-ANG	10





GLASSMAN Anterior Resection Clamp

Product Number	Jaw	Jaw Length (mm)	Length (in)
GS0609	90°,SERR	25	10
GS0610	90°,SERR	50	10
GS0611	90°,SERR	75	10



GS0609



GS0610



GS0611



GS0609 - GS0611

MASTERSON Uterine Clamp

Product Number	Jaw	Length (in)
OG0167	ST	8
OG0168	Cd	8



OG0167



OG0168

OG0167 - OG0168

Z-Type Hysterectomy Clamp



OG0169 - OG0176, OG0181 - OG0184



OG0177 - OG0180, OG0185 - OG0188

Product Number	Jaw	Handle	Length (in)
OG0169	ST	ST	8
OG0170	SLIGHT-Cd	ST	8
OG0171	Cd	ST	8
OG0172	ANG	ST	8
OG0173	ST	ST	9
OG0174	SLIGHT-Cd	ST	9
OG0175	Cd	ST	9
OG0176	ANG	ST	9
OG0177	ST	ANG	9
OG0178	SLIGHT-Cd	ANG	9
OG0179	Cd	ANG	9
OG0180	ANG	ANG	9
OG0181	ST	ST	12
OG0182	SLIGHT-Cd	ST	12
OG0183	Cd	ST	12
OG0184	ANG	ST	12
OG0185	ST	ANG	12
OG0186	SLIGHT-Cd	ANG	12
OG0187	Cd	ANG	12
OG0188	ANG	ANG	12



OG0169, OG0173, OG0177, OG0181, OG0185



OG0172, OG0176, OG0180, OG0184, OG0188



OG0170, OG0174, OG0178, OG0182, OG0186



OG0171, OG0175, OG0179, OG0183, OG0187



Z-Type Hysterectomy Clamp (continued)

Product Number	Jaw	Handle	Length (in)
OG0189	ST	ST	14
OG0190	SLIGHT-Cd	ST	14
OG0191	Cd	ST	14
OG0192	ANG	ST	14



OG0189



OG0190



OG0191




OG0192



OG0189 - OG0192

GELPI Perineal Retractor W/ Ball Stop

Product Number	Length (in)
GS1356	7



GELPI Vaginal Retractor

Product Number	Tip	Length (in)
GS1357	SHARP	5
GS1358	SHARP	6
GS1359	SHARP	7
GS1360	BLUNT	4



RIGBY Vaginal Retractor

Product Number	Length (in)
OG0193	7



DELEE Universal Retractor

Product Number	Length (in)
OG0194	9



OG0194

O' SULLIVAN-O' CONNOR Retractor

Product Number	Description	Blade Size (mm)
OG0195	COMPLETE	185 x 143
OG0196	SMALL BLADE	80 x 149
OG0197	LARGE BLADE	155 x 110

OG0195 Includes:
1 Retractor
2 Blades OG0196 and
1 Blade OG0197



OG0195 - OG0197



OG0196



OG0197



OG0198

COLLIN Vaginal Specula


Product Number	Blade Size (mm)
OG0198	95/85 x 18

CUSCO Vaginal Specula

Product Number	Blade Size (mm)
OG0199	75 x 17



OG0199



OG0200 - OG0202


CUSCO Vaginal Specula

Product Number	Blade Size (mm)
OG0200	80 x 32
OG0201	90 x 36
OG0202	110 x 37



SEYFFERT Vaginal Specula

Product Number	Blade Size (mm)	Length (in)
OG0203	100 x 17	9
OG0204	110 x 33	9



OG0203 - OG0204

SEYFFERT Speculum, 1 Side Open

Product Number	Blade Size (mm)	Length (in)
OG0205	110 x 33	9



OG0205

SEMM Vaginal Specula, 45° Angle

Product Number	Blade Size (mm)	Length(in)
OG0206	90 x 14	7



OG0206

PEDERSON Vaginal Speculum

Product Number	Blade Size (mm)
OG0207	75 x 12.7
OG0208	100 x 22
OG0209	115 x 16
OG0210	120 x 25



OG0207 - OG0210



**PEDERSON Vaginal Speculum, Large,
Insulated Coating, W/ Evacuation Vent**

Product Number	Blade Size (mm)
OG0211	120 x 25



OG0211



OG0212 - OG0215

GRAVES Vaginal Speculum

Product Number	Blade Size (mm)
OG0212	75 x 20
OG0213	95 x 30
OG0214	105 x 33
OG0215	125 x 40
OG0216	102 x 34 W/SIDE OPENING



OG0216

KRISTELLER Vaginal Speculum, Set of 2 each

Product Number	Blade Size (mm)	Length (in)
OG0217	70 x 15	7
OG0218	70 x 27	7
OG0219	75 x 30	7
OG0220	90 x 36	7
OG0221	110 x 20	9
OG0222	110 x 27	9
OG0223	110 x 30	9
OG0224	110 x 36	9

Each blade size comes as a set of 2 - with one upper and lower blade.



OG0217 - OG0224



KRISTELLER Vaginal Specula Upper Blades

Product Number	Blade Size (mm)	Length (in)
OG0225	75 x 10	7
OG0226	75 x 23	7
OG0227	80 x 25	7
OG0228	100 x 32	7
OG0229	115 x 14	8
OG0230	115 x 23	8
OG0231	115 x 26	9
OG0232	120 x 30	9



OG0225 - OG0232



OG0225



OG0226



OG0227



OG0228



OG0229



OG0230



OG0231



OG0232



OG0233 - OG0241

KRISTELLER Vaginal Specula Lower Blades

Product Number	Blade Size (mm)	Length (in)
OG0233	70 x 15	7
OG0234	70 x 27	7
OG0235	75 x 30	7
OG0236	90 x 36	7
OG0237	110 x 18	8
OG0238	110 x 27	8
OG0239	110 x 20	9
OG0240	110 x 30	9
OG0241	110 x 36	9



OG0233



OG0234



OG0235



OG0236



OG0237



OG0238



OG0239



OG0240



OG0241



OG0242

KALLMORGEN Vaginal Specula Set of 2

Product Number	Blade Size (mm)	Length (in)
OG0242	75 x 39	7



OG0243

KALLMORGEN Vaginal Specula Set of 2

Product Number	Blade Size (mm)	Length(in)
OG0243	95 x 39	8



OG0244 - OG0246

SIMS Vaginal Retractor

Product Number	Blade Size (mm)
OG0244	70 x 25
OG0245	76 x 36
OG0246	89 x 38

SIMS Vaginal Speculum, Double Ended

Product Number	Blade Size (mm)	Length (in)
OG0247	69 x 25 / 70 x 30	5
OG0248	70 x 26 / 75 x 30	5
OG0249	80 x 36 / 90 x 40	7



OG0247 - OG0249



JACKSON Vaginal Speculum Lower Blade

Product Number	Blade	Blade Size (mm)	Length(in)
OG0250	SLIGHT-CONCAVE	89 x 38	7
OG0251	SLIGHT-CONCAVE	102 x 38	7



OG0250 - OG0251



OG0250



OG0251

EASTMAN Vaginal Retractor

Product Number	Blade Size (mm)	Length (in)
OG0252	75 x 40	7
OG0253	89 x 38	7
OG0254	100 x 40	7



OG0252 - OG0254

HEANEY Retractor

Product Number	Blade Size (mm)	Length(in)
OG0255	45 x 25	7
OG0256	73 x 29	7
OG0257	83 x 29	7



OG0255 - OG0256

OG0257



OG0258 - Og0260

DOYEN Vaginal Speculum

Product Number	Blade Size (mm)	Length(in)
OG0258	60 x 45	9
OG0259	90 x 45	9
OG0260	112 x 45	9



OG0258



OG0259



OG0260



DOYEN Vaginal Speculum

Product Number	Blade Size (mm)	Length (in)
OG0261	60 x 35	9
OG0262	90 x 35	9
OG0263	120 x 35	9
OG0264	60 x 45	9
OG0265	90 x 45	9
OG0266	120 x 45	9
OG0267	60 x 60	9
OG0268	90 x 60	9
OG0269	120 x 60	9
OG0270	160 x 45	9
OG0271	160 x 60	9



OG0261 - OG0271



OG0272 - OG0283

BREISKY Vaginal Speculum

Product Number	Blade Size (mm)	Length (in)
OG0272	100 x 20	12
OG0273	100 x 25	12
OG0274	100 x 30	12
OG0275	100 x 35	12
OG0276	100 x 40	12
OG0277	130 x 20	12
OG0278	130 x 25	12
OG0279	130 x 30	12
OG0280	130 x 35	12
OG0281	130 x 40	12
OG0282	160 x 40	14
OG0283	180 x 40	14




AUVARD Speculum, W/ Detachable Weight

Product Number	Blade Size (mm)	Length (in)
OG0284	65 x 38-28	9



OG0284



OG0285 - OG0286

AUVARD Speculum, W/ Detachable Weight

Product Number	Blade Size (mm)	Length (in)
OG0285	75 x 42-40	9
OG0286	105 x 43-40	9



Og0287 - Og0289

AUVARD Speculum

Product Number	Blade Size (mm)	Length (in)
OG0287	48 x 70	9
OG0288	85 x 38	10
OG0289	120 x 38	10

STEINER-AUVARD Vaginal Speculum, Malleable, Goose Neck

Product Number	Blade Size (mm)
OG0290	140 x 30



OG0290



DOYEN Lateral Speculum

Product Number	Blade Size (mm)	Length (in)
OG0291	60 x 30	9
OG0292	80 x 30	10

OG0291 - OG0292

SCHERBACK Vaginal Speculum Set

Product Number	Description	Blade Size (mm)
OG0293	SET COMPLETE	-
OG0294	BLADE	80 x 30
OG0295	BLADE	85 x 35
OG0296	BLADE	85 x 40
OG0297	BLADE	55 x 40
OG0298	BLADE	70 x 45
OG0299	HANDLE	-
OG0300	WEIGHT	-



OG0293




OG0294



OG0295



OG0296



OG0297 - OG0298

OG0293 Includes:
1 Handle OG0299
1 Weight OG0300
4 Blades of each
OG0294 - OG0298



TUEBINGER Vaginal Retractor

Product Number	Blade Size (mm)	Length (in)
OG0301	95 x 20	10



OG0301

OG0301

GLENNER Retractor

Product Number	Blade Size (mm)	Length (in)
OG0302	135 x 25	11

OG0302





MARTIN Vaginal Retractor

Product Number	Blade Size (mm)	Length (in)
OG0303	105 x 27	10



SIMON Vaginal Retractor

Product Number	Blade Size (mm)	Length (in)
OG0304	115 x 22	11
OG0305	115 x 27	11



OG0304

OG0305

OG0304 - OG0305



KOGAN Endospeculum

Product Number	Length (in)
OG0306	9



Og0306

Og0306



Og0307

Og0307

KOGAN Cervix/Urethra Endospeculum, W/ Scale, Ratcheted

Product Number	Length (in)
OG0307	9

HEGAR Dilator, Single Ended

Product number	Tip Size	Length (in)
OG0308	1	7
OG0309	2	7
OG0310	3	7
OG0311	4	7
OG0312	5	7
OG0313	6	7
OG0314	7	7
OG0315	8	7
OG0316	9	7
OG0317	10	7
OG0318	11	7
OG0319	12	7
OG0320	13	7
OG0321	14	7
OG0322	15	7
OG0323	16	7
OG0324	17	7
OG0325	18	7
OG0326	19	7
OG0327	20	7
OG0328	21	7



OG0308 - OG0350



HEGAR Dilator, Single Ended (Continued)



OG0308 - OG0350

Product Number	Tip Size (mm)	Length (in)
OG0329	22	7
OG0330	23	7
OG0331	24	7
OG0332	25	7
OG0333	26	7
OG0334	1.5	7
OG0335	2.5	7
OG0336	3.4	7
OG0337	4.5	7
OG0338	5.5	7
OG0339	6.5	7
OG0340	7.5	7
OG0341	8.5	7
OG0342	9.5	7
OG0343	10.5	7
OG0344	11.5	7
OG0345	12.5	7
OG0346	13.5	7
OG0347	14.5	7
OG0348	15.5	7
OG0349	16.5	7
OG0350	17.5	7



OG0351 - OG0352

HEGAR Dilator Set, Single Ended

Product Number	Item
OG0351	SET W/ 14 DILATORS, 4-17mm
OG0352	SET W/ 26 DILATORS, 1-26mm

OG0351 Includes:
Rack for Dilators 4-12 mm
Rack for Dilators 13-17 mm
Metal Case: 200 x 100 x 35 mm

Includes Dilators:
OG0311 - OG0324

OG0352 Includes:
Rack for Dilators 1-17 mm
Rack for Dilators 18-26 mm
Metal Case: 330 x 220 x 80 mm

Includes Dilators:
OG0308 - OG0333



OG0353 - OG0362

HEGAR Uterine Dilator, Double Ended

Product Number	Tip Size (mm)	Length (in)
OG0353	SET OF 8	-
OG0354	1/2	8
OG0355	3/4	8
OG0356	5/6	8
OG0357	7/8	8
OG0358	9/10	8
OG0359	11/12	8
OG0360	13/14	8
OG0361	15/16	8
OG0362	17/18	8

Set Includes:
OG0355 - OG0362



HEGAR Uterine Dilator, Slight Curved, Double Ended

Product Number	Tip Size (mm)	Length (in)
OG0363	SET OF 8	-
OG0364	1/2	8
OG0365	3/4	8
OG0366	5/6	8
OG0367	7/8	8
OG0368	9/10	8
OG0369	11/12	8
OG0370	13/14	8
OG0371	17/18	8
OG0372	17/18	8

OG0363 Set of 8 includes OG0364 - OG0372




OG0363 - OG0372

PRATT Uterine Dilator

Product Number	Tip Size (mm)	Length (in)
OG0373	COMPLETE SET	-
OG0374	13/14	11
OG0375	17/18	11
OG0376	21/23	11
OG0377	25/27	11
OG0378	29/31	11
OG0379	33/35	11
OG0380	37/39	11
OG0381	41/43	11

OG0373 Set of 8 includes OG0374 - OG0381



OG0373 - OG0381

HANK Dilator, French Scale

Product Number	Tip Size (mm)	Length (in)
OG0382	COMPLETE SET	10
OG0383	9/10	10
OG0384	11/12	10
OG0385	13/14	10
OG0386	15/16	10
OG0387	17/18	10
OG0388	19/20	10

OG0382 Set of 6 includes OG0383, OG0384 - OG0388



OG0382 - OG0388



HANK Dilator, American Scale

Product Number	Tip Size (mm)	Length (in)
OG0389	COMPLETE SET	-
OG0390	9/10	11
OG0391	11/12	11
OG0392	13/14	11
OG0393	15/16	11
OG0394	17/18	11
OG0395	19/20	11

OG0389 Set of 6 includes OG0390 - OG0395

OG0389 - OG0395

SIMS Uterine Probe, Graduated (cm)

Product Number	Tip	Tip Size (mm)	Blade	Length (in)
OG0396	RGD	4	Cd	13
OG0397	MALL	4	Cd	13



OG0396 - OG0397

MARTIN Uterine Probe, Graduated (cm)

Product Number	Blade	Length (in)
OG0398	Cd	12



OG0398



OG0399

VALLEIX Uterine Probe, Graduated (cm), W/ Slider

Product Number	Tip Size(mm)	Blade	Length (in)
OG0399	4	Cd	10

SIMPSON Uterine Probe, Graduated (cm)

Product Number	Tip	Tip Size(mm)	Blade	Length (in)
OG0400	MALL	4	Cd	12

OG0400



BRAUN Uterine Depressor, Double Ended

Product Number	Blade Size (mm)	Length (in)
OG0401	24/28	11



OG0401



**ASCH Uterine Secretion Scoop,
Malleable, Double Ended**

Product Number	Length (in)
OG0402	8

**NOVAK Endometrial Biopsy
Suction Curette**

Product Number	Length (in)
OG0403	9



OG0403



**RANDALL Endometrial Suction and
Biopsy Curette**

Product Number	Length (in)
OG0404	9

OG0404

BERLIN Endometrial Biopsy Curette

Product Number	Tip Size (mm)	Length (in)
OG0405	3.5	11



OG0405

KEVORKIAN-YOUNGE Biopsy Curette

Product Number	Tip	Description	Length (in)
OG0406	SLIGHT-Cd	W/O Basket	13
OG0407	SLIGHT-Cd	W/ BASKET	12



OG0406



OG0407



RECAMIER Uterine Curette



Product Number	Tip Size(#)	Tip Size (mm)	Tip	Length (in)
OG0408	00	4.5	SHARP	12
OG0409	0	7.5	SHARP	12
OG0410	1	8.5	SHARP	12
OG0411	2	10.5	SHARP	12
OG0412	3	12.0	SHARP	12
OG0413	4	14.0	SHARP	12
OG0414	5	15.0	SHARP	12
OG0415	6	16.5	SHARP	12
OG0416	7	18.0	SHARP	12
OG0417	8	19.5	SHARP	12
OG0418	9	21.5	SHARP	12
OG0419	10	23.5	SHARP	12
OG0420	1	8.5	SHARP/MALL	12
OG0421	2	10.5	SHARP/MALL	12
OG0422	3	12.0	SHARP/SMALL	12
OG0423	4	14.0	SHARP/MALL	12



OG0408




OG0429




OG0408 - OG0438



OG0410



OG0411




OG0412

RECAMIER Uterine Curette (Continued)


Product Number	Tip Size (#)	Tip Size (mm)	Tip	Length (in)
OG0424	5	15.0	SHARP/SMALL	12
OG0425	6	16.5	SHARP/SMALL	12
OG0426	7	18	SHARP/SMALL	12
OG0427	8	19.5	SHARP/SMALL	12
OG0428	0	7.5	BLUNT	12
OG0429	1	8.5	BLUNT	12
OG0430	2	10.5	BLUNT	12
OG0431	3	12.0	BLUNT	12
OG0432	4	14.0	BLUNT	12
OG0433	5	15.0	BLUNT	12
OG0434	6	16.5	BLUNT	12
OG0435	7	18.0	BLUNT	12
OG0436	8	19.5	BLUNT	12
OG0437	9	21.5	BLUNT	12
OG0438	10	23.5	BLUNT	12
OG0439	1	34.5	BLUNT	14
OG0440	2	39.5	BLUNT	14



OG0439 - OG0440



OG0439



OG0440



SIMS Uterine Curette

Product Number	Tip Size (#)	Tip Size (mm)	Tip	Length (in)
OG0441	1	7.5	SHARP	10
OG0442	2	8.5	SHARP	10
OG0443	3	10.5	SHARP	10
OG0444	4	12.0	SHARP	10
OG0445	5	14.0	SHARP	10
OG0446	6	15.0	SHARP	10
OG0447	1	7.5	SHARP/MALL	10
OG0448	2	8.5	SHARP/MALL	10
OG0449	3	10.5	SHARP/MALL	10
OG0450	4	12.0	SHARP/MALL	10
OG0451	5	14.0	SHARP/MALL	10
OG0452	6	15.0	SHARP/MALL	10
OG0453	1	7.5	SHARP	10

OG0441 - OG0467



1



4




2



5



3



6

SIMS Uterine Curette (Continued)

Product Number	Tip Size (#)	Tip Size (mm)	Tip	Length (in)
OG0454	2	8.5	SHARP	10
OG0455	3	10.5	SHARP	10
OG0456	4	12.0	SHARP	10
OG0457	5	14.0	SHARP	10
OG0458	6	15.0	SHARP	10
OG0459	000	4.5	SHARP/MALL	11
OG0460	00	5.6	SHARP/MALL	11
OG0461	0	6.2	SHARP/MALL	11
OG0462	1	7.3	SHARP/MALL	11
OG0463	2	8.5	SHARP/MALL	11
OG0464	3	9.4	SHARP/MALL	11
OG0465	4	10.6	SHARP/MALL	11
OG0466	5	11.5	SHARP/MALL	11
OG0467	6	12.6	SHARP/MALL	11



OG0441 - OG0467



000



00



0



1



2



3



4



5



6




BUMM Uterine Curette, Blunt Malleable

Product Number	Tip Size (#)	Tip Size (mm)	Length (in)
OG0468	0	16.0	11
OG0469	1	18.0	11
OG0470	2	19.5	11
OG0471	3	21.5	11
OG0472	4	23.5	11
OG0473	5	26.0	11
OG0474	6	28.5	11



0




1



2



3



4



5






6

OG0468 - OG0474

SIMS THOMAS Uterine Curette, Blunt/Malleable

Product Number	Tip Size (#)	Tip Size (mm)	Length (in)
OG0475	0	6.0	11
OG0476	1	8.5	11
OG0477	2	10.5	11
OG0478	3	12.0	11
OG0479	4	14.0	11
OG0480	5	15.0	11
OG0481	6	16.5	11



OG0475 - OG0481

KELLY-GRAY Uterine Curette

Product Number	Tip Size (#)	Tip Size (mm)	Length (in)
OG0482	1	5.0	10
OG0483	2	6.0	10
OG0484	3	7.5	10
OG0485	4	8.5	10



OG0482 - OG0485



OG0482



OG0483



OG0484



OG0485



HEANEY Uterine Curette, Malleable

Product Number	Tip	Length (in)
OG0486	TOOTHED	10



OG0486

OG0486

COLLIN Uterine Curette, Rigid

Product Number	Tip	Length (in)
OG0487	TOOTHED	12



OG0487

BLAKE Uterine Curette, Double Ended

Product Number	Tip	Length (in)
OG0488	SHARP/BLUNT	11




OG0488

BUMM Uterine Curette

Product Number	Tip Size (#)	Tip Size (mm)	Tip	Length (in)
OG0489	1	23.5	SHARP/MALL	12
OG0490	3	28.5	SHARP/MALL	12
OG0491	1	23.5	BLUNT/MALL	12
OG0492	2	26.0	BLUNT/MALL	12
OG0493	3	28.5	BLUNT/MALL	12



OG0489 - OG0493



WALLICH Placenta/Ovum Scoop

Product Number	Tip	Length (in)
OG0494	Cd	17




OG0494

OG0494



Guide Sheath, IOWA Trumpet

Product Number	Tip Size(mm)	Description	Length (in)
OG0495	-	-	7
OG0496	1 x 185	INJECTION NEEDLE, LUER-LOCK HUB	-



Uterine Cannula W/ 2 Cones

Product Number	Length (in)
OG0497	14



Product Number	Page	Product Number	Page	Product Number	Page
OG0001.....	01	OG0026.....	08	OG0053.....	17
OG0002.....	01	OG0027.....	08	OG0054.....	17
OG0003.....	01	OG0028.....	09	OG0055.....	17
OG0004.....	01	OG0029.....	09	OG0056.....	18
OG0005.....	01	OG0030.....	09	OG0057.....	18
OG0006.....	01	OG0031.....	09	OG0058.....	19
OG0007.....	02	OG0032.....	09	OG0059.....	19
OG0008.....	02	OG0033.....	09	OG0060.....	20
OG0009.....	02	OG0034.....	09	OG0061.....	20
OG0010.....	02	GS0282.....	10	OG0062.....	20
OG0011.....	02	GS0350.....	10	OG0063.....	21
OG0012.....	02	GS0351.....	10	OG0064.....	21
OG0013.....	03	OG0035.....	11	OG0065.....	22
OG0014.....	03	OG0036.....	11	OG0066.....	22
OG0015.....	04	OG0037.....	11	OG0067.....	22
OG0016.....	04	OG0038.....	11	OG0068.....	23
OG0017.....	05	OG0039.....	12	OG0069.....	23
OG0018.....	05	OG0040.....	12	OG0070.....	23
OG0019.....	05	OG0041.....	12	OG0071.....	23
OG0020.....	06	OG0042.....	12	OG0072.....	24
OG0021.....	06	OG0043.....	13	OG0073.....	24
OG0022.....	06	OG0044.....	13	OG0074.....	24
OG0023.....	07	OG0045.....	13	OG0075.....	25
OG0024.....	07	OG0046.....	13	OG0076.....	25
GS0214.....	07	OG0047.....	14	OG0077.....	26
GS0176.....	08	OG0048.....	14	OG0078.....	26
GS0177.....	08	OG0049.....	15	OG0079.....	27
GS0178.....	08	OG0050.....	15	OG0080.....	27
GS0179.....	08	OG0051.....	16	OG0081.....	28
OG0025.....	08	OG0052.....	16	OG0082.....	28



Product Number	Page	Product Number	Page	Product Number	Page
OG0083.....	29	OG0113.....	39	OG0143.....	50
OG0084.....	29	OG0114.....	39	OG0144.....	50
OG0085.....	30	OG0115.....	40	OG0145.....	50
OG0086.....	30	OG0116.....	40	OG0146.....	50
OG0087.....	31	OG0117.....	40	OG0147.....	50
OG0088.....	31	OG0118.....	41	OG0148.....	50
OG0089.....	32	OG0119.....	41	OG0149.....	50
OG0090.....	32	OG0120.....	42	OG0150.....	50
OG0091.....	32	OG0121.....	42	OG0151.....	50
OG0092.....	32	OG0122.....	42	OG0152.....	51
OG0093.....	33	OG0123.....	43	OG0153.....	51
OG0094.....	33	OG0124.....	44	OG0154.....	51
OG0095.....	33	OG0125.....	45	OG0155.....	52
OG0096.....	33	OG0126.....	45	OG0156.....	52
OG0097.....	34	OG0127.....	45	OG0157.....	52
OG0098.....	34	OG0128.....	45	OG0158.....	52
OG0099.....	34	OG0129.....	46	OG0159.....	53
OG0100.....	35	OG0130.....	46	OG0160.....	53
OG0101.....	35	OG0131.....	46	OG0161.....	53
OG0102.....	35	OG0132.....	47	OG0162.....	53
OG0103.....	36	OG0133.....	47	OG0163.....	54
OG0104.....	36	OG0134.....	48	OG0164.....	54
OG0105.....	36	OG0135.....	48	OG0165.....	55
OG0106.....	37	OG0136.....	49	OG0166.....	55
OG0107.....	37	OG0137.....	49	GS0609.....	56
OG0108.....	38	OG0138.....	49	GS0610.....	56
OG0109.....	38	OG0139.....	49	GS0611.....	56
OG0110.....	38	OG0140.....	49	OG0167.....	56
OG0111.....	38	OG0141.....	50	OG0168.....	56
OG0112.....	38	OG0142.....	50	OG0169.....	57

Product Number	Page	Product Number	Page	Product Number	Page
OG0170.....	57	OG0195.....	60	OG0225.....	66
OG0171.....	57	OG0196.....	60	OG0226.....	66
OG0172.....	57	OG0197.....	60	OG0227.....	66
OG0173.....	57	OG0198.....	61	OG0228.....	66
OG0174.....	57	OG0199.....	61	OG0229.....	66
OG0175.....	57	OG0200.....	61	OG0230.....	66
OG0176.....	57	OG0201.....	61	OG0231.....	66
OG0177.....	57	OG0202.....	61	OG0232.....	66
OG0178.....	57	OG0203.....	62	OG0233.....	67
OG0179.....	57	OG0204.....	62	OG0234.....	67
OG0180.....	57	OG0205.....	62	OG0235.....	67
OG0181.....	57	OG0206.....	62	OG0236.....	67
OG0182.....	57	OG0207.....	63	OG0237.....	67
OG0183.....	57	OG0208.....	63	OG0238.....	67
OG0184.....	57	OG0209.....	63	OG0239.....	67
OG0185.....	57	OG0210.....	63	OG0240.....	67
OG0186.....	57	OG0211.....	64	OG0241.....	67
OG0187.....	57	OG0212.....	64	OG0242.....	68
OG0188.....	57	OG0213.....	64	OG0243.....	68
OG0189.....	58	OG0214.....	64	OG0244.....	69
OG0190.....	58	OG0215.....	64	OG0245.....	69
OG0191.....	58	OG0216.....	64	OG0246.....	69
OG0192.....	58	OG0217.....	65	OG0247.....	69
GS1356.....	59	OG0218.....	65	OG0248.....	69
GS1357.....	59	OG0219.....	65	OG0249.....	69
GS1358.....	59	OG0220.....	65	OG0250.....	70
GS1359.....	59	OG0221.....	65	OG0251.....	70
GS1360.....	59	OG0222.....	65	OG0252.....	70
OG0193.....	59	OG0223.....	65	OG0253.....	70
OG0194.....	60	OG0224.....	65	OG0254.....	70



Product Number	Page	Product Number	Page	Product Number	Page
OG0255.....	71	OG0285.....	73	OG0315.....	81
OG0256.....	71	OG0286.....	73	OG0316.....	81
OG0257.....	71	OG0287.....	74	OG0317.....	81
OG0258.....	71	OG0288.....	74	OG0318.....	81
OG0259.....	71	OG0289.....	74	OG0319.....	81
OG0260.....	71	OG0290.....	74	OG0320.....	81
OG0261.....	72	OG0291.....	75	OG0321.....	81
OG0262.....	72	OG0292.....	75	OG0322.....	81
OG0263.....	72	OG0293.....	75	OG0323.....	81
OG0264.....	72	OG0294.....	75	OG0324.....	81
OG0265.....	72	OG0295.....	75	OG0325.....	81
OG0266.....	72	OG0296.....	75	OG0326.....	81
OG0267.....	72	OG0297.....	75	OG0327.....	81
OG0268.....	72	OG0298.....	75	OG0328.....	81
OG0269.....	72	OG0299.....	75	OG0329.....	82
OG0270.....	72	OG0300.....	75	OG0330.....	82
OG0271.....	72	OG0301.....	76	OG0331.....	82
OG0272.....	72	OG0302.....	77	OG0332.....	82
OG0273.....	72	OG0303.....	78	OG0333.....	82
OG0274.....	72	OG0304.....	79	OG0334.....	82
OG0275.....	72	OG0305.....	79	OG0335.....	82
OG0276.....	72	OG0306.....	80	OG0336.....	82
OG0277.....	72	OG0307.....	80	OG0337.....	82
OG0278.....	72	OG0308.....	81	OG0338.....	82
OG0279.....	72	OG0309.....	81	OG0339.....	82
OG0280.....	72	OG0310.....	81	OG0340.....	82
OG0281.....	72	OG0311.....	81	OG0341.....	82
OG0282.....	72	OG0312.....	81	OG0342.....	82
OG0283.....	72	OG0313.....	81	OG0343.....	82
OG0284.....	73	OG0314.....	81	OG0344.....	82

Product Number	Page	Product Number	Page	Product Number	Page
OG0345.....	82	OG0375.....	84	OG0416.....	87
OG0346.....	82	OG0376.....	84	OG0417.....	87
OG0347.....	82	OG0377.....	84	OG0418.....	87
OG0348.....	82	OG0378.....	84	OG0419.....	88
OG0349.....	82	OG0379.....	84	OG0420.....	88
OG0350.....	82	OG0380.....	84	OG0421.....	88
OG0351.....	83	OG0381.....	84	OG0422.....	88
OG0352.....	83	OG0382.....	84	OG0423.....	88
OG0353.....	83	OG0383.....	84	OG0424.....	88
OG0354.....	83	OG0394.....	84	OG0425.....	88
OG0355.....	83	OG0395.....	84	OG0426.....	88
OG0356.....	83	OG0396.....	84	OG0427.....	88
OG0357.....	83	OG0397.....	84	OG0428.....	88
OG0358.....	83	OG0398.....	84	OG0429.....	88
OG0359.....	83	OG0399.....	85	OG0430.....	88
OG0360.....	83	OG0400.....	85	OG0431.....	88
OG0361.....	83	OG0401.....	85	OG0432.....	88
OG0362.....	83	OG0402.....	85	OG0433.....	88
OG0363.....	84	OG0403.....	85	OG0434.....	88
OG0364.....	84	OG0404.....	85	OG0435.....	89
OG0365.....	84	OG0405.....	85	OG0436.....	89
OG0366.....	84	OG0406.....	85	OG0437.....	89
OG0367.....	84	OG0407.....	85	OG0438.....	89
OG0368.....	84	OG0408.....	85	OG0439.....	89
OG0369.....	84	OG0409.....	86	OG0450.....	89
OG0370.....	84	OG0410.....	86	OG0452.....	89
OG0371.....	84	OG0411.....	86	OG0453.....	89
OG0372.....	84	OG0413.....	87	OG0454.....	89
OG0373.....	84	OG0414.....	87	OG0455.....	89
OG0374.....	84	OG0415.....	87	OG0456.....	89



Product Number	Page	Product Number	Page	Product Number	Page
OG0436.....	89	OG0466.....	91	OG0496.....	96
OG0437.....	89	OG0467.....	91	OG0497.....	96
OG0438.....	89	OG0468.....	92		
OG0439.....	89	OG0469.....	92		
OG0440.....	89	OG0470.....	92		
OG0441.....	90	OG0471.....	92		
OG0442.....	90	OG0472.....	92		
OG0443.....	90	OG0473.....	92		
OG0444.....	90	OG0474.....	92		
OG0445.....	90	OG0475.....	93		
OG0446.....	90	OG0476.....	93		
OG0447.....	90	OG0477.....	93		
OG0448.....	90	OG0478.....	93		
OG0449.....	90	OG0479.....	93		
OG0450.....	90	OG0480.....	93		
OG0451.....	90	OG0481.....	93		
OG0452.....	90	OG0482.....	93		
OG0453.....	90	OG0483.....	93		
OG0454.....	91	OG0484.....	93		
OG0455.....	91	OG0485.....	93		
OG0456.....	91	OG0486.....	94		
OG0457.....	91	OG0487.....	94		
OG0458.....	91	OG0488.....	94		
OG0459.....	91	OG0489.....	95		
OG0460.....	91	OG0490.....	95		
OG0461.....	91	OG0491.....	95		
OG0462.....	91	OG0492.....	95		
OG0463.....	91	OG0493.....	95		
OG0464.....	91	OG0494.....	95		
OG0465.....	91	OG0495.....	96		

C
A

Product Name	Page
ABEL Tenaculum Forceps.....	19
ADAIR Uterine Forceps.....	14
ASCH Uterine Secretion Scoop.....	87
AUVARD Speculum.....	74
AUVARD Speculum, W/ Detachable Weight.....	73

B

Product Name	Page
BARNES-NEVILLE Obstetrical Forceps.....	35
BARRETT Placenta Forceps.....	38
BARRETT Uterine Tenaculum Forceps.....	15
BERGER Biopsy Forceps.....	46
BERLIN Endometrial Biopsy Curette.....	85
BIERER Ovum Forceps.....	38
BILL Traction Handle.....	35
BLAKE Uterine Curette.....	94
BOZEMANN Forceps.....	13
BOZEMANN Uterine Scissors.....	4
BOZEMANN-DOUGLAS Forceps.....	14
BRAUN Tenaculum Forceps.....	16
BRAUN Uterine Depressor.....	86
BRAUN-STADLER Episiotomy Scissors.....	6
BREISKY Vaginal Speculum.....	72
BUMM Uterine Curette.....	92, 95

Product Name
Page

COLLIN Uterine Curette.....	28
COLLIN Uterine Elevating Forceps.....	28
COLLIN Vaginal Specula.....	61
CUSCO Vaginal Specula.....	61
CZERNY Uterine Tenaculum Forceps.....	20

D
Product Name
Page

DARTIGUES Elevating Forceps.....	29
DE LEE Obstetrical Forceps.....	33
DELEE Universal Retractor.....	60
DOYEN Lateral Speculum.....	75
DOYEN Myoma Screw.....	31
DOYEN Ovum Forceps.....	27
DOYEN Vaginal Speculum.....	71, 72
DUPLAY Tenaculum Forceps.....	18

E
Product Name
Page

EASTMAN Vaginal Retractor.....	70
ELLIOTT Obstetrical Forceps.....	33
EPPENDORF Biopsy Forceps.....	48

F
Product Name
Page

FLETCHER Sponge Forceps.....	25
FLETCHER-VAN DOREN Uterine Polyp, Sponge Holding Forceps.....	24



G

Product Name	Page
GELPI Perineal Retractor W/ Ball Stop.....	59
GELPI Vaginal Retractor.....	59
GLASSMAN Anterior Resection Clamp.....	56
GLENNER Retractor.....	77
GORDON Uterine Vulsellum Forceps.....	20
GRAVES Vaginal Speculum.....	64
GREEN-ARMYTAGE Forceps.....	37
Guide Sheath, IOWA Trumpet.....	96

H

Product Name	Page
HANK Dilator, American Scale.....	85
HANK Dilator, French Scale.....	84
HEANEY Hysterectomy Forceps.....	50, 51, 52
HEANEY Hysterectomy Forceps, Heavy.....	51
HEANEY Needle Holder.....	10
HEANEY Retractor.....	71
HEANEY Uterine Curette, Malleable.....	94
HEANEY-BALLENTINE Hysterectomy Forceps.....	50
HEANEY-REZEK Hysterectomy Forceps.....	49
HEGAR Dilator Set, Single Ended.....	83
HEGAR Dilator, Single Ended.....	81
HEGAR Uterine Dilator, Double Ended.....	83
HENROTIN Uterine Vulsellum Forceps.....	19
HEYWOOD-SMITH Polyp Forceps.....	26
HOLZBACH Hysterectomy Forceps.....	49
HULKA-KENNWICK Uterine Elevating Forceps.....	30

I

Product Name	Page
IOWA Membrane Puncturing Forceps.....	30

J

Product Name	Page
JACKSON Vaginal Speculum Lower Blade.....	70
JACOBS Vulsellum Forceps.....	22
JORGENSEN Scissors.....	3

K

Product Name	Page
KALLMORGEN Vaginal Specula Set of 2.....	68
KELLY-GRAY Uterine Curette.....	93
KEVORKIAN-YOUNGE Biopsy Curette.....	87
KEVORKIAN-YOUNGE Biopsy Forceps.....	47
KIELLAND Obstetrical Forceps.....	32
KOGAN Cervix/Urethra Endospeculum.....	80
KOGAN Endospeculum.....	80
KRISTELLER Vaginal Specula Lower Blades.....	67
KRISTELLER Vaginal Specula Upper Blades.....	66

L

Product Name	Page
LITTLEWOOD Tenaculum Forceps.....	22
LUIKART-KIELLAND Obstetrical Forceps.....	34
LUIKART-MCLANE Obstetrical Forceps.....	34
LUIKART-SIMPSON Obstetrical Forceps.....	34
LUIKHART Obstetrical Forceps.....	36

M

Product Name	Page
MAIER Forceps.....	11
MARTIN Uterine Probe.....	85
MARTIN Vaginal Retractor.....	78
MASTERSON Uterine Clamp.....	56
MAYO Dissecting Scissors.....	8
MCCLINTOCK Placenta and Ovum Forceps.....	43
MCLANE-TUCKER Obstetrical Forceps.....	36
METZENBAUM Wavecut Dissecting Scissors.....	7
MURLESS Head Extractor.....	36
MUSEUX Tenaculum Forceps.....	22

N

Product Name	Page
NAEGELE Obstetrical Forceps.....	32
NOTO Polyp Forceps.....	26
NOVAK Endometrial Biopsy Suction Curette.....	87

O

Product Name	Page
O' SULLIVAN-O' CONNOR Retractor.....	60
Obstetrical Forceps.....	32
OVERSTREET Endometrial Polyp Forceps.....	25

P

Product Name	Page
PEDERSON Vaginal Speculum.....	63
PELKMANN Forceps.....	12
PHANEUF Hysterectomy Forceps.....	49
PIPER Obstetrical Forceps, F/ Breach.....	33
POZZI Tenaculum Forceps.....	15, 18
PRATT Tenaculum Forceps.....	20, 21
PRATT Uterine Dilator.....	84

R

Product Name	Page
RANDALL Endometrial Suction and Biopsy Curette...	87
RECAMIER Uterine Curette.....	88, 89
RIGBY Vaginal Retractor.....	59

S

Product Name	Page
SAENGER Placenta and Ovum Forceps.....	42
SCHERBACK Vaginal Speculum Set.....	75
SCHROEDER Tenaculum Forceps.....	16
SCHUHMACHER Umbilical Cord Scissors.....	6
SCHUMACHER Biopsy Forceps, S-Shape.....	45
SCHUMACHER Biopsy Rongeur.....	45
SELLHEIM Obstetric Lever.....	37
SEMM Vaginal Specula, 45° Angle.....	62
SEYFFERT Speculum, 1 Side Open.....	62
SEYFFERT Vaginal Specula.....	62
SIEBOLD Uterine Scissors.....	3
SIMON Vaginal Retractor.....	79
SIMPSON C-Section Forceps.....	35
SIMPSON Uterine Probe.....	86
SIMS Scissors.....	1
SIMS THOMAS Uterine Curette.....	93
SIMS Uterine Curette.....	90, 91
SIMS Uterine Probe.....	85
SIMS Uterine Scissors.....	8
SIMS Uterine Wavecut Scissors.....	7



S

Product Name	Page
SIMS Vaginal Retractor.....	69
SIMS Vaginal Speculum.....	69
SIMS-SIEBOLD Uterine Scissors.....	4
SOMER Uterine Elevating Forceps.....	29
STEINER-AUVARD Vaginal Speculum.....	74

T

Product Name	Page
TEALE Tenaculum Forceps.....	21
THOMS-GAYLOR Biopsy Forceps, S-Shape.....	44
TISCHLER Biopsy Alligator Forceps.....	47
TISCHLER Cervical Biopsy Punch Forceps.....	46
TUEBINGER Vaginal Retractor.....	76

U

Product Name	Page
Umbilical Cord Scissors.....	5
Uterine Cannula W/ 2 Cones.....	96
Uterine Scissors.....	9

V

Product Name	Page
VALLEIX Uterine Probe.....	86

W

Product Name	Page
WALLICH Placenta/Ovum Scoop.....	95
WERTHEIM Atraumatic Clamp.....	55
WERTHEIM Atraumatic Parametrium Clamp.....	53
WERTHEIM Hysterectomy Forceps.....	54
WERTHEIM Parametrium Clamp.....	54
WERTHEIM Scissors.....	2, 8
WERTHEIM-CULLEN Clamp.....	55
WIENER Hysterectomy Forceps.....	52
WINTER Placenta and Ovum Forceps.....	39, 40, 41
WITTNER Biopsy Forceps.....	48

Z

Product Name	Page
Z-Type Hysterectomy Clamp.....	57, 58



Scan to know more

Thomson Surgicals Pvt. Ltd

31/26, 1st Floor, 6th Cross, Bhavani Nagar,
Bengaluru - 560 029, Karnataka, India

P +91 80 26685076 M +91 94808 12041

E mail@thomsonsurgicals.com

Mfg Units Bengaluru - Karnataka | Palakad - Kerala

W www.thomsonsurgicals.com

